NORDIC BIOLABS PRESENTERAR

# Transfusionsnytt

1/2019

NYHET!

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
DaraEx inhibits anti-CD38 ("Dara"). Erythrocytes were preincubated with DaraEx or Tris buffered saline (TBS) as incubation control. Subsequently, the cells were incubated with anti-CD38 ("Dara") (left), anti-CD38 + anti-Kell serum (middle), anti-CD38 + anti-Fy(a) serum (right) and analyzed in indirect antiglobulin test (IAT).

Dara treated cells react false positive in the absence of antibodies; this effect is diminished after DaraEx treatment (left). Dara treated cells react false positive in the presence of anti-Kell (middle) or anti-Fy(a) serum (right); this effect is diminished after DaraEx treatment.

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# Daratumumab Interference in Red Cell Antibody Detection: Specific Inhibition with DaraEx

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## Introduction

The anti-cancer drug Daratumumab, an anti-CD38 antibody, has recently been approved by the European Commission for treatment of multiple myeloma. CD38 is a surface molecule expressed on a number of tissues and cells, including erythrocytes.

The high titers of Daratumumab in patient sera lead to a strong interference in indirect antiglobulin tests (IAT), where most, if not all, of the reactions turn positive. This can happen up to 6 month after the last treatment with Daratumumab. The workarounds that have been tested so far include

- the removal of CD38 on the erythrocytes by treating them with DTT or Trypsin,
- inhibition of Daratumumab with soluble CD38 or an anti-idiotypic antibody, and
- matching by genotyping.

All of these approaches have considerable drawbacks and technical difficulties. We therefore decided to test a reagent, called “DaraEx”, which promises to specifically inhibit the reaction of Daratumumab with erythrocytes.

## Materials and Methods

DaraEx is an experimental reagent obtained from imusyn GmbH & Co. KG, Hanover, Germany. Figure 1 shows the protocol for DaraEx treatment that has been used in this study.

## Results


In order to evaluate the efficiency and specificity of DaraEx, we tested a number of cells with Daratumumab, sera, or mixtures of Daratumumab and sera. The reaction of Daratumumab with the erythrocytes was completely abolished when the cells were pre-treated with DaraEx. “Hidden” antibodies (here: anti-s, anti-Kell and anti-Fy(a)) could be revealed, indicating that DaraEx is specific for Daratumumab only (Figure 2 and 3).

As Daratumumab titers in patients can be different from the concentrations used in the previous experiments, we tried to inhibit the reaction of test cells with sera from a patient treated with Daratumumab. Apart from a “mixed-field” result (a line of erythrocytes directly on top of the gel bed), which most likely is not due to Daratumumab, the reaction of the patient sample with the test erythrocytes was inhibited (Figure 4).

## Conclusions

The number of patients treated with Daratumumab will increase strongly in the foreseeable future, bringing new challenges to immunohematology laboratories. The drawbacks of the current workarounds make the search for alternatives very attractive.

Our results show that DaraEx can specifically inhibit Daratumumab induced agglutination of erythrocytes. Other antibodies such as anti-Kell are not affected, making DaraEx more specific than the commonly used DTT-treatment. Further studies are required to fully evaluate the potential of DaraEx with patient sera


**Fig. 1: Protocol for DaraEx Treatment.**

Flow chart of the DaraEx protocol. After adding 0.3 vol DaraEx to the erythrocyte solution, the cells were incubated for 30 min at room temperature on a laboratory shaker. Subsequently, the cells were used in the IAT.


**1 imusyn GmbH & Co. KG, Feodor-Lynen-Str. 5, 30625 Hannover, Germany**

**2 Institute for Transfusion Medicine, Hannover Medical School, Carl-Neuberg-Str. 1, 30625 Hannover, Germany**

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
**Fig. 2: Specific Inhibition of Daratumumab by DaraEx.** Cell no. 7 (MNS-blood group: M, N, s) and cell no. 11 (N, S) of the ID-DiaPanel (Bio-Rad) were pre-incubated either with TBS or DaraEx as indicated. The erythrocytes were then incubated with either 300 µg/ml Daratumumab (DARA) in 0.9% NaCl, an anti-s serum (anti-s), or an anti-s serum spiked with 300 µg/ml Daratumumab (DARA + anti-s) in an IAT. Shown are the results of an ID-Gel Card (Bio-Rad).


**Fig. 3: KEL and FY are not Affected by DaraEx.** Cell no. 2 (KEL-blood group: K, k, FY-blood group: Fy(b)) and 4 (k, Fy(a)) of the ID-DiaPanel were pre-incubated either with TBS or DaraEx as indicated. The erythrocytes were then incubated with either an anti-K serum (anti-Kell), an anti-Fy(a) serum (anti-Fy(a)), an anti-K serum spiked with 300 µg/ml Daratumumab (DARA + anti-Kell), or an anti-Fy(a) serum spiked with 300 µg/ml Daratumumab (DARA + anti-Fy(a)) in an IAT. Shown are the results from ID-Gel Cards.

**Fig. 4: Inhibition with DaraEx in a Patient Sample.** Cells no. 9 and 11 of the ID-DiaPanel were pre-incubated either with DaraEx or TBS as indicated. The erythrocytes were then incubated with either 300 µg/ml Daratumumab (DARA) in 0.9% NaCl or a serum of a patient treated with Daratumumab in an IAT. Shown are the results from an ID-Gel Card.

#### Acknowledgements

We thank Hein Hustinx of the BSD SRK Bern for kindly providing Daratumumab patient sera.

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Nordic Biolabs AB, Box 7293, 187 14 Täby,