TPP Techno Plastic Products AG
Zollstrasse 7, CH-8219 Trasadingen, Switzerland
Tel. +41 (0)52 687 01 87, Fax +41 (0)52 687 01 77
info@tpp.ch, www.tpp.ch

Gebrauchsanleitung TPP - Zellkulturflaschen / *Instructions for Use TPP Tissue Culture Flasks*

Anleitung

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Instructions

TPP tissue culture flasks are intended for cell culture in manual use. The growth area is opto-mechanically treated for an optimal cell adhesion and proliferation. The chanted neck decreases the risk of media contact at the inside of the screw cap. It also enables a better access for pipettes and scrapers to reach the corners of the flask.

- *Open product packaging and inspect for intact bag and product.*
- *Fill medium: Introduce medium or inoculate cells with a pipette alongside the wall, or direct injection on the growth area of the flask.*
- *Close flask:*

1. *VENT screw cap for open culture systems, by opening the cap the gas exchange takes place*
 - o *Introduce screw cap until it clicks in the VENT position.*
 - ✓ *Position for defined aeration of medium: rectangle on 12 o'clock.*
 - *Further turn screw cap over the click about 90° to close flask gas tight.*
 - ✓ *Position for gas tight: rectangle on 3 o'clock*
2. *Filter screw cap with integrated hydrophobic membrane of pore size 0.22 µm for sterile gas exchange after the cap is screwed down tightly.*
 - o *Screw cap down tightly.*

- *Incubate: Stapling ring and rim at the angled neck result in a stable and secure stand in the incubator.*
- *Harvest cells:*
 - o *Harvest adherent cells by trypsinisation or scape with cell scraper.*
 - o *Harvest suspension cells by exchange of medium.*

TPP Techno Plastic Products AG
Zollstrasse 7, CH-8219 Trasadingen, Switzerland
Tel. +41 (0)52 687 01 87, Fax +41 (0)52 687 01 77
info@tpp.ch, www.tpp.ch

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Advice and Recommendation
For homogenous cultivation of the cells it is most important to level the incubator and the continuous supply of medium. Sources of vibrations must be eliminated to avoid circular growth pattern. The necessary gas exchange takes place thought the screw cap.

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Technical Data

Materialien

Flasche	PS
Kappe	PE
Membrane	PTFE, Porengrösse 0.22 µm

Materials

Flask	*PS*
Screw cap	*PE*
Membrane	*PTFE, pore size 0.22 µm*

Abmessungen	***Measurement***	**90025**	**90026**	**90075**	**90076**	**90150**	**90151**	**90300**	**90301**
Kappe	*Cap*	VENT		VENT		VENT		VENT	
Kappe	*Cap*		Filter		Filter		Filter		Filter
Membrane µm	*Membrane µm*		0.22		0.22		0.22		0.22
Höhe mm	*Heights mm*	29		40		50		50	
Breite mm	*Width mm*	51		87		122		140	
Länge mm	*Length mm*	92		155		210		275	
Wachstums-fläche cm^2	*Growth area cm^2*	25		75		150		300	
Empf Volumen ml	*Rec. volume ml*	3 – 5		8 – 15		15 – 30		30 – 50	
Max Volumen ml	*Max. volume ml*	15		65		165		410	