

NORDIC BIOLABS PRESENTS THE

Christmas Calendar 2018

The offers in the Christmas Calendar are valid only one specific day and can not be combined with other offers or discounts from Nordic Biolabs. Please refer to the Christmas Calendar in order to get the right price when ordering. VAT and shipping cost will be added, but no small order fee or mosa will be charged.

Order tel. 08 630 85 10, e-mail info@nordicbiolabs.se

More information about the products can be found in our webshop - www.nordicbiolabs.se.

**3
dec**

Flask-trap aspirator, volume 1L

FTA-1 Flask-trap aspirator is designed for aspiration/removal of alcohol/buffer remaining quantities from microtest tube walls during DNA/RNA purification and other macromolecule reprecipitation techniques.

FTA-1, Flask-trap aspirator
Art no: BS-040108-AAG

Today! 3 065 SEK
List price 4 026 SEK



**4
dec**

Sprout® Mini Centrifuge is compact, easy to use and comes with a host of extras

- Small footprint and quiet operation
- Palm-shaped lid
- Change out rotors with the push of a button

Maximum Speed: approx. 6,000 rpm
Maximum RCF: 2000 x g

First come ...

Art no	Description	List price SEK	Now! SEK
HS120301	Sprout Mini-Centrifuge,Green	3 069	1 535
HS120396	Sprout Mini-Centrifuge,Purple	3 069	1 535
HS120395	Sprout Mini-Centrifuge,Yellow	3 069	1 535



Healthrow Scientific
Design & Innovation

**5
dec**

Linear Stand for a wide range of pipette brands

When not in use, the pipette should be stored in an upright position in order to avoid contamination from desks and benches. Linear Stand is designed for all Sartorius mechanical and electronic pipettes and is also compatible with a wide range of other pipette makes.

Linear Stand for mechanical
and electronic pipettes
Art no: 725620

Today! 510 SEK/pcs
List price 849 SEK/pcs



sartorius

**6
dec**

Vortex V-1 plus - the ultimate choice for mixing

Two modes:

1. Continuous operation;
2. Impulse operation
(activated by pressing the cap with the tube's bottom)

V-1 Personal vortex
Art no: BS-010203-AAG

Today! 1 325 SEK
List price 1 899 SEK



bioSan
Medical, Biological
Research & Technology

 **nordic biolabs**

**7
dec**

Adapt-A-Rack - flexible Multi-Tube Rack



- 5 mL–50 mL tubes sizes
- autoclavable
- 2-pack



Art no	Description	List price SEK	Now! SEK
HS120184	Tube Rack, Combo Blue/Green	111	60
HS120189	Multi-Tube Rack, Pink	111	60
HS120187	Tube Rack, Blue	111	60
HS120185	Tube Rack, Combo Blue/Yellow	111	60

**10
dec**

TruCool® Ice pans and buckets



Made of ethylene-vinyl acetate (EVA) foam, great non-toxic alternative to conventional ice pans. Highly insulative, stackable, unbreakable, lightweight and leak-proof. Ideal for use with ice, dry ice (-78°C), liquid nitrogen (-196°C), alcohol or saline solutions, warm solutions (+93°C).

First come ...



Art no	Description	List price SEK	Now! SEK
BCS-111PK	Ice Pan, maxi 9L, pink	2 016	1 318
BCS-111PL	Ice Pan, maxi 9L, purple	2 016	1 318
BCS-113GR	Ice Pan, midi 4L, lime green	1 695	992
BCS-113PK	Ice Pan, midi 4L, pink	1 470	992
BCS-113PL	Ice Pan, midi 4L, purple	1 470	992
BCS-115-25G	Ice Bucket, round, 2,5L, green	1 470	745
BCS-115-25PK	Ice Bucket, round, 2,5L, lime green	1 092	745
BCS-115GR	Ice Bucket, round 4L, lime green	1 300	844
BCS-115PK	Ice Bucket, round 4L, pink	1 198	844
BCS-115PL	Ice Bucket, round 4L, purple	1 198	844

**11
dec**

Tacta mechanical pipette with superb comfort and reliability

The forces needed for a pipetting cycle, from tip attachment to pipetting and finally tip ejection, are exceptionally low with Tacta, protecting you from WRULD (Work Related Upper Limb Disorder).



Art no		List price SEK	Now! SEK
Tacta Mechanical Pipette, 1-ch			
LH-729010	0.1-3 ul	2 812	1 658
LH-729020	0.5-10 ul	2 812	1 658
LH-729030	2-20 ul	2 812	1 658
LH-729050	10-100 ul	2 812	1 658
LH-729060	20-200 ul	2 812	1 658
LH-729070	100-1000 ul	2 812	1 658

**12
dec**

Picus electronic pipettes - highly sophisticated and completely ergonomic design

Highest Level of Ergonomics and reliable results. Uniquely lightweight (100 g) and compact design for effortless pipetting experience. Conveniently located, soft-touch operating button and electronic tip ejection help minimize muscle strain.



Art no	Description	List price SEK	Now! SEK
Picus Electronic Pipette, 1-ch			
735021	0.2-10 ul, AC-adaptor Universal	8 531	5840
735041	5-120 ul, AC-adaptor Universal	8 531	5840
735061	10-300 ul AC-adaptor Universal	8 531	5840
735081	50-1000ul, AC-adaptor Universal	8 531	5840
Picus Electronic Pipette, 8-ch			
735321	0.2-10 ul, AC-adaptor Universal	13 256	8995
735361	10-300 ul, AC-adaptor Universal	13 256	8995
735391	50-1200 ul AC-adaptorUnivers.	13 256	8995
Picus Electronic Pipette, 12-ch			
735441	5-120 ul, AC-adaptor Universal	15 772	10760
735461	10-300 ul, AC-adaptor Universal	15 772	10760
735491	50-1200 ul, AC-adaptorUniversal	15 772	10760

**13
dec**

Programmable Rotator Multi Bio RS-24 performs several motion types in one module

Programmable Rotator Multi Bio RS-24 performs several motion types in one module. Microprocessor control allows performing not only Vertical rotation (1) of the platform, but also Reciprocal rotation (2) as well as Vibration (3).

Multi Bio RS-24
Art no: BS-010117-AAG

Today 4 167 SEK
List price 5 556 SEK



**14
dec**

CoolCell- freeze your cells without alcohol



Ensure standardized controlled-rate -1°C/minute cell freezing in a -80°C freezer - without alcohol or any fluids. Proven for use with a variety of cell types including stem cells, primary cells and cell lines. Easy to use - simply fill with cryogenic vials and place in a -80°C freezer. For 1.0 mL - 5.0 mL cryogenic vials, and injectable serum vials for cell therapy applications.

biocision

Art no	Description	List price SEK	Now! SEK
BCS-405	CoolCell LX, purple, 12 tubes	2 280	1 465
BCS-405G	CoolCell LX, green, 12 tubes	2 280	1 465
BCS-405PK	CoolCell LX, pink, 12 tubes	2 280	1 465
BCS-170	CoolCell FTS30, purple, 30 tubes	5 998	3 975
BCS-170PK	CoolCell FTS30, pink, 30 tubes	5 998	3 975
BCS-170G	CoolCell FTS, green, 30 tubes	5 998	3 975
BCS-170O	CoolCell FTS, orange, 30 tubes	5 998	3 975

**17
dec**

Spark free refrigerators and freezers for blood banks, drugs/ vaccines and/or laboratories

- Spark free interior
- Refrigerator: Automatic defrosting
- Freeze: Manual defrosting
- Lockable door
- Digital Display with max/min temp-log function
- Light and sound alarm at high/ low temp and open door



LABCOLD
LABORATORY APPLIANCE

Art.nr.	Beskrivning	List price SEK	Now! SEK
RLPR0214	Spark free Refrigerator 49L, bench model, Lock	7 250	5 275
RLPR0517	Spark free Refrigerator 150L, undercounter, Lock	8 810	6 410
RLPR1514	Spark free Refrigerator 400L, Free standing	19 500	14 200
RLVF0214	Spark free Freezer 40L, bench model, Lock	8 740	6 360
RLVF1514	Apark free Freezer 400L, Free standing	21 600	15 725

**18
dec**

TruCool® Hinged Cryoboxes - Hinged lid, water-resistant fiberboard



Lid is easy to open when frozen and stays attached to base, minimizing contamination risk. Recommended for refrigerators, freezers, or ultra-low freezer storage. Store up to 81 tubes. Dimensions (L x W x H): 13.3 x 13.3 x 5.1 cm.

biocision

Art no	Description	List price SEK	Now! SEK
BCS-206B	Hinged CryoBox, blue 5pk	522	265
BCS-206G	Hinged CryoBox, green 5pk	522	265
BCS-206MC	Hinged CryoBox, multi, 5pk	522	265
BCS-206O	Hinged CryoBox, orange 5pk	522	265
BCS-206P	Hinged CryoBox, purple 5pk	522	265
BCS-206PK	Hinged CryoBox, pink 5pk	522	265

19
dec

The strongest pipette boy on the market



43%
discount

- **Versatil usage:** fits all commercially available glass and plastic serological pipettes in the volume range 1-100 ml
- **Speed control:** an adjustment thumb wheel for setting the pipetting speed
- **Gravity dispensing:** dispensing drop-by-drop by applying a slight pressure on the dispense button

Pipettor Turbo Fix
Art no: 94700

Today! 2 275 SEK

List price 3 990 SEK



20
dec

Easy dispensing with easy-moving dosing button - Ripette® pro

- **Flexible:** With the Ripette® pro, you can work flexibly and practically with 12 tip sizes, 10 volume adjustments and 120 program steps.
- **Clear:** Volume adjustments and dosing volumes for different tips are shown in a vision panel on the Ripette pro.
- **Resistance:** By using selected premium materials the Ripette pro shows a high resistance towards a broad range of chemicals.

Ripette pro
Art no: 42020-000

Today! 2 000 SEK

List price 3 688 SEK

46%
discount



ritter

21
dec

Combi-Spin, Mini-Centrifuge/Vortex

Mini-centrifuge/vortex Combi-Spin FVL-2400N is specially designed for genetic engineering research (for PCR-diagnostics experiments).

Combi-Spin provides simultaneous mixing and separation of samples, using centrifuge and mixing modules, located on the common spin-module.



34%
discount

FVL-2400N Combi-Spin, Mini-Centrifug/Vortex
Art no: BS-010202-AAA

Today 2 395 SEK

List price 3 630 SEK



Merry Christmas
and A Happy New Year
From all of us at Nordic Biolabs

All offers are valid one specific day and can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please refer to the Christmas Calendar in order to get the right price when ordering. With reservation for misprints.

 **nordic biolabs**
www.nordicbiolabs.se

Nordic Biolabs AB, Box 7293, 187 14 Täby,
Tel vxl 08-630 85 00, Fax 08-630 09 05, info@nordicbiolabs.se