

NORDIC BIOLABS PRESENTS THE

# Christmas Calendar 2017

**The offers in the Christmas Calendar are valid only one specific day** and can not be combined with other offers or discounts from Nordic Biolabs. Please refer to the Christmas Calendar in order to get the right price when ordering. VAT and shipping cost will be added, but no small order fee or mosa will be charged.

Order tel. 08 630 85 10, e-mail [info@nordicbiolabs.se](mailto:info@nordicbiolabs.se)

More information about the products can be found in our webshop - [www.nordicbiolabs.se](http://www.nordicbiolabs.se).

**1  
dec**

## PBS for cell culture work

Buffer cell culture with GE Healthcare HyClone Dulbecco's Phosphate Buffered Saline (PBS), solution.

**30%  
discount**



GE Healthcare  
Life Sciences

Art no	Description	List price SEK	Now! SEK
SH30028.FS	D-PBS (1X) without Ca, Mg, 500mLx6	426	<b>299</b>
SH30378.02	D-PBS 10X without Ca, Mg, 1L	184	<b>130</b>

**4  
dec**

## Sprout® Mini Centrifuge is compact, easy to use and comes with a host of extras

- Small footprint and quiet operation
- Palm-shaped lid
- Change out rotors with the push of a button

**Maximum Speed:** approx. 6,000 rpm

**Maximum RCF:** 2000 x g



Healthrow Scientific  
Design & Innovation

Art no	Description	List price SEK	Now! SEK
HS120301	Sprout Mini-Centrifuge, Green	3 069	<b>2 000</b>
HS120396	Sprout Mini-Centrifuge, Purple	3 069	<b>2 000</b>
HS120395	Sprout Mini-Centrifuge, Yellow	3 069	<b>2 000</b>

**5  
dec**

## Linear Stand for a wide range of pipette brands

When not in use, the pipette should be stored in an upright position in order to avoid contamination from desks and benches. Linear Stand is designed for all Sartorius mechanical and electronic pipettes and is also compatible with a wide range of other pipette makes.

Linear Stand for mechanical and electronic pipettes  
Art no: 725620

**Today! 500 SEK/pcs**

List price 849 SEK/pcs



**41%  
discount**

sartorius

**6  
dec**

## Vortex V-1 plus - the ultimate choice for mixing

*Two modes:*

1. Continuous operation;
2. Impulse operation (activated by pressing the cap with the tube's bottom)

V-1 Personal vortex

Art no: BS-010203-AAG

**Today! 1 300 SEK**

List price 1 809 SEK



**28%  
discount**

bioSan  
Medical, Biological  
Research & Technology

 **nordic biolabs**

**7  
dec**

**TS-100C, Thermo-Shaker with cooling for microtubes and PCR plates**



**35%  
discount**

Thermo-Shaker TS-100C provides intensive mixing and temperature control. Feature examples - fast reaching of specified mixing speed and maintenance of equal amplitude of rotation throughout the Thermo-Shaker block, stability of maintaining the preset temperature in a wide range throughout the Thermo-Shaker's block surface, LCD display indicates preset and current values of temperature, speed and time of operation. Functions of heating and mixing can be performed both simultaneously and independently.

TS-100C, Thermo-Shaker with cooling for microtubes and PCR plates  
Art no: BS-010143-AAI

**Today 6 715 SEK**  
*List price 10 331 SEK*



**8  
dec**

**TruCool® Ice pans and buckets**



**30%  
discount**

Made of ethylene-vinyl acetate (EVA) foam, these pans and buckets are a great non-toxic alternative to conventional ice pans. The containers are highly insulative, stackable, unbreakable, lightweight and leak-proof. Ideal for use with ice, dry ice (-78°C), liquid nitrogen (-196°C), alcohol or saline solutions, or warm solutions (+93°C).

First come ...



Art no	Description	List price SEK	Now! SEK
BCS-111PL	Ice pan, maxi 9L, purple	2 016	1 411
BCS-113PK	Ice pan, midi 4L, pink	1 470	1 029
BCS-113PL	Ice pan, midi 4L, purple	1 470	1 029
BCS-115PK	Ice bucket, round 4L, pink	1 198	838
BCS-211B	Ice pan, mini, 1L, blue	1 028	720

**11  
dec**

**Compact and userfriendly High-Speed Multi Plate Shaker**

**37%  
discount**



High-Speed Multi Plate Shaker MPS-1 can be used in virtually any application by providing adjustable mixing of reagents in microtest plates, PCR plates, deepwell plates and test tubes (shaking tubes 0.2 to 2 ml and vortexing any volume up to 50 ml). MPS-1 features a head for vortexing a single tube.

Shaker can be used in cold rooms or incubators, operating at ambient temperature range +4°C to +40°C. Low voltage external power supply (12 V) provides electrical safety in humid environment.

MPS1 High Speed Multi Plate Shaker  
Art no: BS-010216-A11

**Today 5 700 SEK**  
*List price 9 115 SEK*

**12  
dec**

**Up to 40% discount on our most popular pipettes**

**Tacta**

The forces needed for a pipetting cycle, from tip attachment to pipetting and finally tip ejection, are exceptionally low with Tacta, protecting you from WRULD (Work Related Upper Limb Disorder).

**Picus**

Highest Level of Ergonomics and reliable results. Uniquely lightweight (100 g) and compact design for effortless pipetting experience. Conveniently located, soft-touch operating button and electronic tip ejection help minimize muscle strain.



Picus Tacta



Pipette	List price SEK	Now! SEK
Tacta, mechanical pipette, single channel	2 812	1 799
Tacta, mechanical pipette, 8 channel	7 587	4 855
Tacta, mechanical pipette, 12 channel	8 467	5 419
Picus, electronic pipette, single channel	8 531	5 400
Picus, electronic pipette, 8 channel	13 256	7 900
Picus, electronic pipette, 12 channel	15 772	9 500

**13  
dec**

## Programmable Rotator Multi Bio RS-24 performs several motion types in one module

Microprocessor control allows performing not only Vertical rotation of the platform, but also Reciprocal rotation as well as Vibration. These three motion types can be performed separately, pairwise and in cycles, periodically repeating the sequence of three motion types. Multi-Rotation option of Biosan instruments substantially expands possibilities and enhances efficiency of sample preparation for the examined materials and allows adjusting the mixing procedure according to the individual tasks.

Multi Bio RS-24  
Art no: BS-010117-AAG

**Today 3 200 SEK**  
List price 5 291 SEK



**15  
dec**

## Spark free refrigerators and freezers for blood banks, drugs/ vaccines and/or laboratories

- Spark free interior
- Refrigerator: Automatic defrosting
- Freeze: Manual defrosting
- Lockable door
- Digital Display with max/min temp-log function
- Light and sound alarm at high/ low temp and open door



**LABCOLD**  
LABORATORY REFRIGERATORS & FREEZERS

Art.nr.	Beskrivning	List price SEK	Now! SEK
RLPR0214	Spark free Refrigerator 49L, bench model	9 628	5 500
RLPR0517	Spark free Refrigerator 150L, undercounter	11 736	7 000
RLPR1514	Spark free Refrigerator 400L, Free standing	25 454	15 000
RLVF0214	Spark free Freezer 40L, bench model	14 800	8 500
RLVF0416	Spark free Freezer 124L, undercounter	17 300	10 000
RLVF1514	Apark free Freezer 400L, Free standing	33 000	19 000

**14  
dec**

## CoolCell- freeze your cells without alcohol



Ensure standardized controlled-rate -1°C/minute cell freezing in a -80°C freezer - without alcohol or any fluids. Proven for use with a variety of cell types including stem cells, primary cells and cell lines. Easy to use - simply fill with cryogenic vials and place in a -80°C freezer. For 1.0 mL - 5.0 mL cryogenic vials, and injectable serum vials for cell therapy applications.

**bio**cision

Art no	Description	List price SEK	Now! SEK
BCS-405	CoolCell LX, purple, 12 tubes	2 280	1 700
BCS-405G	CoolCell LX, green, 12 tubes	2 280	1 700
BCS-405PK	CoolCell LX, pink, 12 tubes	2 280	1 700
BCS-406	CoolCell 5 ml LX, 12 tubes	3 190	2 300
BCS-170	CoolCell FTS30, purple, 30 tubes	5 998	4 200
BCS-170PK	CoolCell FTS30, pink, 30 tubes	5 998	4 200
BCS-170G	CoolCell FTS, green, 30 tubes	5 998	4 200
BCS-170O	CoolCell FTS, orange, 30 tubes	5 998	4 200

**18  
dec**

## TruCool® Hinged Cryoboxes - Hinged lid, water-resistant fiberboard



Lid is easy to open when frozen and stays attached to base, minimizing contamination risk. Recommended for refrigerators, freezers, or ultra-low freezer storage. Store up to 81 tubes.

Dimensions (L x W x H): 13.3 x 13.3 x 5.1 cm.

**bio**cision

Art no	Description	Size	List price SEK	Now! SEK
BCS-206O	Hinged CryoBox, orange 5pk	5	522	365
BCS-206PK	Hinged CryoBox, pink 5pk	5	522	365
BCS-206MC	Hinged CryoBox, multi, 5pk	5	522	365
BCS-221G	Hinged CryoBox, LN drain holes, Green	50	4 863	3 404
BCS-221P	Hinged CryoBox, LN drain holes, Purple	50	4 863	3 404
BCS-215P	Hinged CryoBox, 3,5" 81-place, LN drain holes, purple	5	775	542



19  
dec

## The strongest pipette boy on the market



32%  
discount

- **Versatil usage:** fits all commercially available glass and plastic serological pipettes in the volume range 1-100 ml
- **Speed control:** an adjustment thumb wheel for setting the pipetting speed
- **Gravity dispensing:** dispensing drop-by-drop by applying a slight pressure on the dispense button

Pipettor Turbo Fix  
Art no: 94700

**Today! 2 200 SEK**

List price 3 231 SEK



20  
dec

## Easy dispensing with easy-moving dosing button - Ripette® pro

- **Flexible:** With the Ripette® pro, you can work flexibly and practically with 12 tip sizes, 10 volume adjustments and 120 program steps.
- **Clear:** Volume adjustments and dosing volumes for different tips are shown in a vision panel on the Ripette pro.
- **Resistance:** By using selected premium materials the Ripette pro shows a high resistance towards a broad range of chemicals.

Ripette pro  
Art no: 42020-000

**Today! 2 000 SEK**

List price 3 688 SEK



46%  
discount



21  
dec

## Combi-Spin, Mini-Centrifuge/Vortex

Mini-centrifuge/vortex Combi-Spin FVL-2400N is specially designed for genetic engineering research (for PCR-diagnostics experiments).

Combi-Spin provides simultaneous mixing and separation of samples, using centrifuge and mixing modules, located on the common spin-module.



38%  
discount

FVL-2400N Combi-Spin, Mini-Centrifuge/Vortex  
Art no: BS-010202-AAA

**Today 2 150 SEK**

List price 3 457 SEK



22  
dec

## Adapt-A-Rack - flexible Multi-Tube Rack



46%  
discount

- 5 mL–50 mL tubes sizes
- autoclavable
- 2-pack



Art no	Description	List price SEK	Now! SEK
HS120184	Tube Rack, Combo Blue/Green	111	60
HS120189	Multi-Tube Rack, Pink	111	60
HS120187	Tube Rack, Blue	111	60
HS120185	Tube Rack, Combo Blue/Yellow	111	60
HS120190	Multi-Tube Rack, Yellow	111	60



Merry Christmas and A Happy New Year  
From all of us at Nordic Biolabs

All offers are valid one specific day and can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please refer to the Christmas Calendar in order to get the right price when ordering. With reservation for misprints.

 **nordic biolabs**  
[www.nordicbiolabs.se](http://www.nordicbiolabs.se)

Nordic Biolabs AB, Box 7293, 187 14 Täby,  
Tel vx1 08-630 85 00, Fax 08-630 09 05, [info@nordicbiolabs.se](mailto:info@nordicbiolabs.se)