



PharmaMedia Dr. Müller GmbH
Gustav-Throm-Str. 1
D-69181 Leimen
Phone: +49 (0) 62 24-92 80 7-0
Email: info@pmm-leimen.de
www.pmm-leimen.de

Certificate of Origin

PMM Swab

Reference: 885.0120

Lot No: 008821

Production date: 2024-Dec-02

Preliminary expiry date: 2025-Jul-25

Storage: 2 – 25°C

Typical composition per L:

Casein Peptone	17,00 g	Soy Peptone	3,00 g
NaCl	5,00 g	K ₂ HPO ₄	2,50 g
Glucose-D(+) x H ₂ O	2,50 g	Polysorbate 80	5,00 g

This medium can be adjusted and / or supplemented according to the performance criteria required.

Raw Material	Specifications	Results	
Casein Peptone	Tissue	milk	conform
	Animal Source	bovine	conform
	Country of origin	NZ, USA	conform
	Infectivity Category	IB	conform
	Tissue	pancreas	conform
	Animal Source	porcine	conform
	Country of origin	USA	conform
	Infectivity Category	n.a.	conform

n.a.: Not applicable

Infectivity Category according to note for guidance on minimizing the risk of transmitting animal spongiform encephalopathy via human and veterinary medicinal products (EMA/410/01 rev.3)

All other raw materials used for the manufacturing of this product are not of animal origin. The origin of all raw materials used by PharmaMedia Dr. Müller GmbH is certified by our suppliers.

Date: 2025-Feb-03


Signed: Anja Rueder
Quality Control