NORDIC BIOLABS PRESENTS THE

Christmas Calendar 2022

All offers in the Christmas Calendar are valid December 1st to 22nd .

Please think of the environment and order all your products at one occasion.

The offers can not be combined with other offers or discounts from Nordic Biolabs. Please refer to the Christmas Calendar for correct price when ordering. VAT and shipping cost will be added, but no small order fee or mosa will be charged.

Order tel. 08 630 85 10, e-mail info@nordicbiolabs.se

More information about the products can be found in our webshop - www.nordicbiolabs.se.

1 dec

Vortex V-1 plus - the ultimate choice for mixing

30% discount



Two modes:

1. Continuous operation
2. Impulse operation (activated by pressing the cap with the tube's bottom)

V-1 Personal vortex
Art no: BS-010203-AAG

Now! 1750 SEK
List price 2531 SEK

biosan Medical–Biological Research & Technologies

2 dec

Nordic Biolabs Mini Centrifuge is compact, easy to use and comes with a host of extras

50% discount

- Small footprint and quiet operation
- Palm-shaped lid
- Change rotors with the push of a button

Maximum Speed: approx. 6,000 rpm
Maximum RCF: 2000 x g



Sprout Plus Mini-Centrifuge, Blue
Art no: HS120611

Now! 2050 SEK
List price 4158 SEK

Heathrow Scientific design | function | innovation

5 dec

FTA-1, Flask-trap aspirator, volume 1L

Designed for aspiration/removal of residual liquids (buffers, alcohols etc) from cell culture vessels, test tubes or microtubes, in cell culture or molecular biology (ex DNA/RNA precipitation).



30% discount

FTA-1, Flask-trap aspirator Volume 1L
Art no: BS-040108-AAG

Now! 3690 SEK
List price 5272 SEK

biosan Medical–Biological Research & Technologies

6 dec

PDS-250 is ready-to-use for eliminating DNA and RNA from suitable surfaces

27% discount

The use of PDS-250 both before and after PCR analysis is fast, easy and ideal to maintain a clean work area and thereby saves time and expenses.



PDS-250
Art no: BS-040107-DK

Now! 375 SEK
List price 515 SEK

biosan Medical–Biological Research & Technologies

nordic biolabs

www.nordicbiolabs.se info@nordicbiolabs.se ordertel 08-630 85 10 fax 08-630 09 05

7 dec

Restock with high quality water and buffer! Your processes deserve it, not only during Christmas

30% discount




cytiva

CAPRICORN SCIENTIFIC

Art no.	Description	List price SEK	Now! SEK
Water For Injection, WFI			
SH31191.02	LM-HyPure WFI Quality Water	88	62
SH31191.03	LM-HyPure WFI Quality Water	126	88
Dulbecco's Phosphate Buffered Saline, D-PBS			
PBS-1A	D-PBS (1x) w/o Ca Mg w/o Phenol Red 500 mL	115	81
PBS-10XA	D-PBS (10x) w/o Ca Mg w/o Phenol Red 500 mL	226	158
PBS-2A	D-PBS (1x) w Ca Mg w/o Phenol Red 500 mL	112	78

8 dec

Sterile universal filter tips in rack



50% discount

Universal filter tips, RNase-, DNase-, Metal- and Pyrogen free. Good compatibility with major pipette brands like Eppendorf, Gilson etc.

- Sterilized by gamma irradiation SAL10-6 (ISO11137)
- Manufactured in a Class 100,000 cleanroom environment under ISO 13485:2016 and ISO 9001:2015 quality management system

deltalab

Art no.	Description	List price SEK	Now! SEK
D151020	20 µl Filter tip, in rack, 96x20 tips	2290	1145
D151100	100 µl Filter tip, in rack, 96x20 tips	2290	1145
D153200	200 µl Filter tip, in rack, 96x20 tips	2290	1145
D151000	1000 µl Filter tip, in rack, 96x20 tips	2781	1391

9 dec

Oakton pH 550 Benchtop pH Meter Kit with Probe, Stand, and pH Buffers

Economical meter for benchtop pH measurement.

- Automatically recognizes calibration solutions
- Compact design saves valuable benchtop space
- Data logging allows you to store and recall up to 50 groups of data
- IP54-rated housing is dust- and spill-proof
- Bundled kit for convenience and savings

Oakton pH 550 Benchtop pH Meter Kit
with Probe, Stand, and pH Buffers
Art no: EW-35419-28

Now! 7130 SEK
List price 8900 SEK

20% discount



Cole-Parmer
an antylia scientific company

12 dec

The strongest pipette boy on the market



30% discount

- **Versatil usage:** fits all commercially available glass and plastic serological pipettes in the volume range 1-100 ml
- **Speed control:** an adjustment thumb wheel for setting the pipetting speed
- **Gravity dispensing:** dispensing drop-by-drop by applying a slight pressure on the dispense button

Pipettor Turbo Fix
Art no: 94700

Now! 3600 SEK
List price 5170 SEK

TPP

13 dec — Centrifuge tubes with benefits

- White writing surface on tube cap, works with any writing instrument
- Caps that cannot be over-tightened
- Graduations all the way down to 0.5 ml



- Continuous thread, external moulded graduations in 0,5 ml increments. Frosted area (55 x 10 mm).
- Blue cap in polyethylene with hermetical closure thanks to its internal elastic obturation
- Not autoclavable.
- Centrifugations force: 5000 xg



TPP

deltalab

Art no.	Description	List price SEK	Now! SEK
91050	Centrifuge tube, 50 ml, conical bottom, PP, 30x115 mm, 10x20 pcs	1627	1100
429920-DL	Centrifuge tube, 15 ml, conical bottom, PP, sterile, 5x100 pcs	1140	790

14 dec — 3D "Sunflower Mini-Shaker" with adjustable three-dimensional smooth rotation of the platform

26% discount



The mini-shaker is designed for mixing blood samples, for minigel staining and destaining, sample washing, blot hybridization reactions.

3D, Sunflower Mini-Shaker
Art no: BS-010151-AAG

Now! 3295 SEK
List price 4454 SEK

biosan
Medical–Biological Research & Technologies

15 dec — Stands for pipettes to avoid contamination from benches and tables



Linear Stand for mechanical and electronical pipettes. The pipette stands hold your pipettes in an upright position while not in use to avoid contamination from benches and tables.

Charging Stand for Picus and eLINE 1 Place Pipette with Universal AC-Adaptor.

SARTORIUS

Art no.	Description	List price SEK	Now! SEK
725620	Linear stand for pipettes	929	604
730981	Charging stand for electronic pipette	3137	2195

16 dec — Portable spill kit that comes in a convenient dispenser box



The dispenser box contains 10 pads (41 x 51 cm). Designed for small spills of chemicals and hazardous liquids. Does not degrade or cause a dangerous chemical reaction with the absorbed liquid. Reusable lightweight plastic box that can be refilled once the sorbents have been used. Perfect to clean up chemical spills or for use in laboratories.

BRADY

Art no.	Size	List price SEK	Now! SEK
813910	10 pads/Kit	755	650
813842	100 pads	1877	1599

19 dec

Riplate PCR plate - affordable high quality

Riplate PCR plate i tested and works well on the widely used CFX96 DX System.

- 96 wells V-bottom, 150µl/well, skirted
- Material: PP (wells) / PC (frame)
- Standard format for use in automation
- Slightly raised well rims for better sealing
- Alphanumeric markings

30% discount



Riplate PCR plate 96, clear thin wall, low profile, 50 pcs/box
Art no: 43001-0305

Now! 2100 SEK
List price 3000 SEK

ritter MEDICAL

20 dec

Extract DNA in 3-8 min

The QuickExtract™ DNA Extraction Solution can be used to rapidly and efficiently extract PCR-ready genomic DNA from almost any sample type using a simple, one-tube protocol that takes only 3-8 minutes.

The QuickExtract method allows for the inexpensive processing of one to hundreds of samples simultaneously, without centrifugation, spin columns or the use of any toxic organic solvent. The method is also compatible with robotic automation.



QuickExtract DNA Extraction Solution 1.0
Art no: QE09050

Now! 4322 SEK
List price 5763 SEK

BIOSEARCH TECHNOLOGIES
Lucigen

21 dec

Spill control products for the lab



1. Kit for small spills of chemicals and hazardous liquids.
2. SpillFix granular is completely natural, producing no carcinogens or dust and has more than four times the absorption power of clay granular.
3. A chemical sorbent soc, filled with chemical meltblown polypropylene that absorbs chemicals and hazardous liquids.

BRADY

Art no.	Description	List price SEK	Now! SEK
813868	Chemical disposable spill kit	1187	998
149265	SpillFix Granular 3 Ltr. Jar	1144	995
813850	Soc, Ø 7.6 cm x 122 cm	1475	1199

22 dec

Reagent reservoirs, 25 and 50 ml



35% discount

- V-shaped well design – max reagent recovery
- Graduation marking with fill line on inner wall
- Drip free pour-off spouts in the corners for spill free pouring

Heathrow Scientific design | function | innovation

Art no.	Material	Sterile	Volume	Pack size	List price SEK	Now! SEK
HS120640	PS	Yes	25 ml	5x20 pcs	2140	1390
HS120642	PS	No	25 ml	100 pcs	861	560
HS120643	PS	Yes	50 ml	5x40 pcs	2516	1640
HS120645	PS	No	50 ml	100 pcs	951	620
HS120646	PVS	Yes	50 ml	5x20 pcs	1704	1105

We wish you a Merry Christmas and a Happy New Year

The offers in the Christmas Calendar are valid until December 23. Please think of the environment and order all your products at one occasion. The offers can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please refer to the Christmas Calendar in order for correct price when ordering. With reservation for misprints.

nordic biolabs
www.nordicbiolabs.se
Nordic Biolabs AB, Box 7293, 187 14 Täby,
Tel vxl 08-630 85 00, Order 08-630 85 10, info@nordicbiolabs.se