NORDIC BIOLABS AB PRESENTS

Products for Life Science

Promo code: YearEnd24

4/2024

Make the nitrile glove choice easy using colour coding!



ONE RISK, ONE COLOUR, ONE GLOVE

- SHIELDskin CHEM™ gloves: protection against CHEMICALS in the laboratory, meeting the ISO 374-1:2016+A1:2018 requirements and the highest expectation of users for comfort & dexterity.
- SHIELDskin™ gloves: protection against BIOHAZARDS and low chemical exposure in the laboratory and meeting the ISO 16604:2004 requirements regarding barrier resistance to virus penetration.
- ecoSHIELD™ gloves: an ECOLOGICAL & ECONOMIC protection in the laboratory for the best compromise between personal and environmental protection.

Select your gloves, adapted to the risk you face!
By choosing the right protection based on risk you will in the end also save the environment and money.

Art no	Description	Length	Thickness	Price/pack when ordering >3 boxes, SEK
67623X	SHIELDskin orange nitrile 260, 10x90 st/fp	260 mm	0,14	1 660
67625X	SHIELDskin orange nitrile 300, 10x50 st/fp	300 mm	0,14	1 365

SHIELD Scientific
compliance comfort protection

X = 1 (XS), 2 (S), 3 (M), 4 (L)

Replace X with the number that corresponds to the size you wish to order.

nordic biolabs

www.nordicbiolabs.se info@nordicbiolabs.se ordertel 08-630 85 10

LIA high speed microcentrifuge

Refrigerated or ventilated small centrifuge with speed rates up to 15 000 rpm and with 21 380 RCF accommodating a high number of samples.

More rotors are available for 1.5/2.0 mL tubes, 5.0 mL tubes and PCR-strips.
Pre-cooling program ensures optimal protection of the sample before warming, ideal for temperature sensitive samples.

Specifications:	
Max Speed	15 000 rpm
Max RCF	21 380 g
Max Loading	18 x 5 mL
Dimensions (HxWxD)	V: 35 x 44 x 45 cm / R: 23 x 33 x 55 cm
Weight	V: 12 kg / R: 30 kg



LIA Centrifuge Ventilated with angular rotor 24x2 mL

Art no: AFI-MICROV-E

Year End Price! 21 600 SEK



LIA Centrifuge Refrigerated with angular rotor 24x2 mL

Art no: AFI-MICROR-E

Year End Price! 56 300 SEK



NuAire CO_2 incubator

Save 33% on all incubators.

- Choose between 160L or 200 L
- Different setups like O_2 and /or humidity control to meet most applications
- ISO Class 5 Clean room environment
- IR sensor for the best environment control
- Lots of options; for instance copper interior and split inner door

Example:

Art no	Description	List price SEK	Right now! SEK
NU5710E-SSRG-N	160L CO2 Incubator with Decon Cycle	116 500	78 050

All offers, except the Promega offer, are valid until December 31, 2024, and can not be combined with other offers or agreements with Nordic Biolabs AB.
All prices are ex VAT. Please refer to promo code "YearEnd24" in order to get the right price when ordering.



STRETCH YOUR YEAR-END BUDGET!

Save 25% with Promegas 4f3 Promo

Promegas popular "Buy 3, Get 1 Free" promotion is back!
Stretch your year-end budget further and save 25% on Promega reagents.

Hurry, these offers are only valid until Monday 16th December.
This promotion is valid on all Promegas standard catalog products, including Maxwell kits, FuGENE products, instruments, and Glo-products. Please note that it cannot be combined with other promotions or special agreements.

How to receive your free product?

Step 1: Combine your orders with at least 3 of the same Promega product
Step 2: Put promo code "4f3" on your order.
Step 3: Please ensure your order is received latest by Monday, December 16, 2024.
Step 4: We'll send you 4 products for the price of 3, saving you 25% on your year-end budget.
Plus, you can easily repeat the process if you need more products!

In other words: *Order 3 identical Promega products and get a 4th one for free! Want more? Order 6 and we'll send you 8! Just make sure your order includes promo code 4F3, and the more you order, the more you save—while supplies last.*



Thermo-shaker made to shake and heat 1 to 4 microtiter plates

- The standard version can be set for heating up to 60°C and is available in 2 designs - one for 2 plates and one for 4 plates.
- The PST-100HL is designed to maintain temperatures up to 100°C.
- Can be operated as 3 independent devices: Incubator, Microplate shaker, Thermo–Shaker
- Speed: 250-1200 rpm



PST-60HL
Shaking and thermostating
2 standard 96-well microplates
Heating <60°C
Outer dimensions: 270x260x125 mm

Art no: BS-010119-AAI
Year End Price! 13 600 SEK
List price: 19 380 SEK



PST-60HL-4
Shaking and thermostating
4 standard 96-well microplates
Heating <60°C
Outer dimensions: 380x390x140 mm

Art no: BS-010128-AAI
Year End Price! 17 800 SEK
List price: 25 340 SEK



PST-100HL
Shaking and thermostating
2 standard 96-well microplates
Heating: 25-100°C
Outer dimensions: 270x260x125 mm

Art no: BS-010142-AAI
Year End Price! 16 200 SEK
List price: 23 110 SEK

Plate Washer

For washing various types of standard 96-well microtiter plates, microstrips as well as microarrays on FastFRAME (rectangular well shape).

Fully programmable ensuring multi-step solution ripening, aspiration (aspiration, combination of aspiration/liquid dispensing and soaking, as well as soaking cycle during a particular period of time). Dispense system of liquid dosage for each channel separately.



- Dispensing volume: 25 ul - 1600 ul, with increments of 25 ul
- 50 programs, nb of washing cycles: 1-15
- Time of one plate wash (300 ul): 45 s

Includes: 3 x 1 liter solution bottles, 1 x 2 liter waste bottle, 1 aerosol collection bottle

Art no: BS-060102-AAI
Year End Price! 40 800 SEK
List price: 58 130 SEK

biosan
Medical–Biological Research & Technologies

All offers, except the Promega offer, are valid until December 31, 2024, and can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please refer to promo code "YearEnd24" in order to get the right price when ordering.

MOLECULAR DEVICES



SPECIAL OFFER!

SpectraMax® Mini Multi-Mode Microplate Reader:

Advanced microplate analysis for academic research

Boost your research with the affordable SpectraMax Mini Multi-Mode Microplate Reader!

Key features

- **Multi-Mode Detection:** Absorbance, fluorescence, and luminescence measurements in one compact instrument.
- **Temperature Control:** Consistent assay outcomes ensured by a built-in system that maintains a stable environment.
- **PathCheck Sensor Technology:** Automatic correction for pathlength variations, delivering reliable, accurate results.
- **SoftMax® Pro Data Acquisition & Analysis Software:** Streamlined data analysis and reporting for efficient research.

Benefits

High-performance

Advanced optics and detection technologies provide precision required for your sensitive assays.

Space-saving

The compact design ensures easy integration into academic lab settings, optimizing available space.


Versatile

Suitable for a variety of applications, including drug discovery, protein quantification, and enzyme kinetics, supporting diverse research goals.

Description	Year End Price SEK
SpectraMax Mini with absorbance, luminescense & fluorescence	167 100
SpectraMax Mini with absorbance & luminescence	157 500
SpectraMax Mini with absorbance & fluorescence	157 500

Plant-based 96-well plates

The 96-well plates are made 100 % from renewable biopolymer and have a 50% lower carbon footprint than that of standard polystyrene plates.



Available in black, white or transparent with a flat bottom. The plates are labeled, non-sterile, without lid and made from polylactic acid (PLA)

- Temperature resistant between -80°C and +50°C
- Centrifugation stability up to 4000 x g
- Working volume 50 – 200 µL/well

Green Elephant BIOTECH



FBS – Choose what is best for your cells

By testing the FBS and reserving the number of bottles you need you can find the best option for your cells.

We offer 2-3 different batches and reserve a minimum of 20 bottles of each batch. Define which works best for your cells. By having both Cytiva and Capricorn as suppliers and also different versions of FBS you have a bigger chance of finding one that works best with you specific cells.



FBS Standard

- South America, USA
- Triple 0.1 µm filtered
- Virus- & bacteria-tested
- Test multiple batches at once
- Batch reservation while testing
- Common cell line-approved

- Natural FBS
- Consitent quality
- Low endotoxin
- Works with all cell types

+ Traditional FBS: Suitable for all cell types

− Pricing is highly dependant on FBS market fluctuations



FBS Xtra

- South America, USA
- Triple 0.1 µm filtered
- Virus- & bacteria-tested
- Test multiple batches at once
- Batch reservation while testing
- Common cell line-approved

- Reduced serum content
- Chemically defined additives
- Reduced lot-to-lot variation
- More sustainable ingredients

+ Best for tumor cells & fast-growing cell types

− Evaluation recommended for long-term cell growth



FBS Advanced

- South America
- Triple 0.1 µm filtered
- Virus- & bacteria-tested
- Test multiple batches at once
- Batch reservation while testing
- Common cell line-approved

- Cost-efficient FBS
- Natural, defined additives
- Minimal lot-to-lot variation
- No further batch testing

+ Best for tumor cells & less demanding cell types

− Transfection optimization may be required



FBS Minis

- South America, USA
- Triple 0.1 µm filtered
- Virus- & bacteria-tested
- Test multiple batches at once
- Batch reservation while testing
- Common cell line-approved

- Pre-aliquoted FBS (10 x 50 ml)
- All serum types are available
- Less risk of contamination
- Reduced plastic waste

+ Less Work. More Time. Easy inventory management

− Slightly higher pricing, but saves valueable time and resources!

Apart from untreated FBS products, we offer most of our Sera with various treatment options, such as heat inactivation, gamma irradiation, exosome depletion, and many others.

All offers, except the Promega offer, are valid until December 31, 2024, and can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please refer to promo code "YearEnd24" in order to get the right price when ordering.

Cell Culture Reagents at year end prices

Classical media, buffers and supplements for cell culture work.

Buy 10 or more from the list below (mix as you like) and get our special year end price!

Product	Volume	Cytiva Art no.	Year End Price Cytiva, SEK	Capricorn Art no.	Year End Price Capricorn, SEK
Buffer and Salt Solutions					
Dulbecco's PBS (10x), w/o Ca & Mg, w/o Phenol Red	500 ml	SH30378.02	141	PBS-10XA	187
Dulbecco's PBS (1x), w/o Ca & Mg, w/o Phenol Red	500 ml	SH30028.01	40	PBS-1A	84
Classical Media					
DMEM High Glucose (4.5 g/l), with L-Glutamine	500 ml	SH30022.01	87	DMEM-HA	87
DMEM High Glucose (4.5 g/l), with L-Glutamine, with Sodium Pyruvate	500 ml	SH30243.01	87	DMEM-HPA	87
RPMI 1640, with L-Glutamine	500 ml	SH30027.01	87	RPMI-A	87
RPMI 1640, w/o L-Glutamine	500 ml	SH30096.01	87	RPMI-XA	87
Growth Supplements					
L-Glutamine (200 mM)	100 ml	SH30034.01	87	GLN-B	87
Antibiotics					
Penicillin/Streptomycin (100x)	100 ml	SV30010	102	PS-B	102
Detachment					
Accutase Cell Detachment Solution*	100 ml	-		ACC-1B	700
Trypsin-EDTA (0.05 %) in DPBS (1x)	100 ml	SH30236.01	79	TRY-1B	95

**The advantages of Accutase over the traditional Trypsin/EDTA treatment are that it is less damaging to cells, thus leading to increased viability, and no inactivation of the enzymes by inhibitors or serum is required.*

cytiva CAPRICORN SCIENTIFIC

SimSorb - a mat for collecting spilled liquids

The mat is a must for work in all laminar flow benches.
No more need for wipes to absorb and remove spills.

- Highly absorbant with a 3-Layer construction
- Latex free
- Non-slip bottom layer
- Low-linting
- Provides a defence against leakages
- Recommendation: ISO 5
- GMP A
- Sterilized by gamma irridation



PROTECT 2CLEAN

Art no.	Description	Mat size	Sterile	Pack Size	List price SEK	Year End Price! SEK
370458	Chemo Prep Mat SimSorb Sterile	280 x 430 mm	Sterile	50 pcs	1 240	960
370459	Chemo Prep Mat SimSorb Sterile	558 x 406 mm	Sterile	50 pcs	2 115	1 509


No contamination using DNA/RNA bench-top UV-cleaner box

No extra lab for your sensitive DNA samples? Problem with contaminated samples?

The DNA/RNA bench-top UV-cleaner box is designed for safe handling of DNA samples.
In this way you only need a small bench space instead of a whole lab to protect your samples from contamination.

Choose between different sizes, materials and see prices below.

biosan
Medical–Biological Research & Technologies

Specifications:				
				
Model	UVC/T-AR (compact)	UVC/T-M-AR (compact)	UVT-B-AR (compact)	UVT-S-AR (double)
Art no	BS-040102-AAA	BS-040104-A06	BS-040109-A06	BS-040107-AAA
List price	*29 600 SEK*	*37 350 SEK*	*31 710 SEK*	*52 850 SEK*
Year End price	**21 400 SEK**	**26 900 SEK**	**22 900 SEK**	**38 100 SEK**
Walls material	Plexiglas: Polymethyl methacrylate (ALTUGLAS EX)	Rear: stainless steel Sides and front: glass (EUROGLASS, Germany)	Rear: stainless steel, Sides: steel with chemical resistant powder coating Front: glass (EUROGLASS, Germany)	Rear: stainless steel Sides and front: glass (EUROGLASS, Germany)
Working surface material	Steel with chemicals resistant powder coating	Stainless steel		
Open UV-lamp	1 × 25W built in bactericidal			2 × 30W built in bactericidal
UV radiation level	15 mW / cm2 / sec			
Radiation type	UV (λ = 253.7 nm), ozone-free			
Digital time setting of direct UV exposure	1 min–24 hrs / non–stop (increment 1 min)			
UV–recirculator	1 × 25 W (efficiency >99% per 1 hour)			1 × 30 W (efficiency >99% per 1 hour)
Daylight lamp	1 × TLD-15W			1 × TLD-30W
Working area dim	645 × 490 mm			1210 × 500 mm
Safety features	Automatic open UV-lamp switch off when screen is open			
Power outlets inside the unit	Inlet for power cords	Inlet for power cords and1 Built-in socket, max. 1,000 W		3 Built-in sockets max. 1,000 W and inlet for power cords
Overall dim W×D×H	720 × 535 × 555 mm	700 x 580 x 555 mm		1250 × 600 × 590 mm

All offers, except the Promega offer, are valid until December 31, 2024, and can not be combined with other offers or agreements with Nordic Biolabs AB.
All prices are ex VAT. Please refer to promo code "YearEnd24" in order to get the right price when ordering.



Pipette Trade-In

*Upgrade your pipetting and save up to 40% on your new pipettes**

The campaign applies to Tacta manual pipettes, single packs

- Ergonomic handle for superior comfort
- Low force when pipetting and tip-ejecting
- Large, clear, user-facing volume display
- Easy to clean and fully autoclavable

Recycle

- We will properly recycle or dispose your old pipettes, in accordance with local regulations
- We will collect and recycle your old pipettes for free

** Based on listprice*

How It Works



1. Select Old Pipettes for Trade-in
Any pipette brand or model will be accepted.

2. Order Trade-in Pipettes
Order Tacta Pipettes on our website - www.nordicbiolabs.se.
Please note number of pipettes you are trading in on the check out page (in the message field) and we will get back to you with information on how to proceed.


3. Free Return Service for Old Pipettes
Old pipettes must be returned within 30 days after receiving new pipettes.



Are your pipettes in need of service and/or calibration?

Let the Pipette Doctor at Nordic Bioabs calibrate your pipettes for optimal accuracy and reliability. Good service, quality-assured calibrations, and repairs included, except for spare parts. We also offer accredited calibration according to ISO 17025.

Please contact us for more information - 08 630 85 05, pipett@nordicbiolabs.se
or visit our website - www.nordicbiolabs.se / Pipettdoktorn.



NGS SERVICES

RNA-Seq, DNA-Seq and bioinformatic solutions

End-to-end or start at any step from extraction to data analysis and reporting, to fit your individual needs.

Why choose to work with Lexogen´s service provider?

Consultation and sharing expertise and support from start to finish.
Customization since specialized in non-standard projects and adapt workflows to your specific needs.
Challenging Material is no problem providing excellent data even from the trickiest of samples.
Consistency and Reproducibility by offering spike-in RNA controls to monitor and evaluate RNA-Seq experiments and data quality.
Data Protection is an important point to consider but since Lexogen works from Austria they operate under the highest data protection regulations.

A selection of what can be offered:

- Gene expression Profiling
- Whole Transcriptome Sequencing
- FFPE RNA Sequencing
- Single-cell RNA sequencing
- Ultra-low input RNA sequencing
- DNA Sequencing

LEXOGEN

Sterile Alcohol 70% High Quality wipes

The wipes are suitable for use in all cleanroom areas as required.

- Each wipe is 28 x 38 cm in size allowing a large surface area to be covered.
- Optimised fluid dosing in conjunction with the highly absorbent 69g per m^2 polyester/cellulose blend enables the product to hold sufficient fluid for in-use requirements without drying out too quickly.
- Packed in thick PET laminate with a resealable label.



Sterile Alcohol (DE) Wipes - 10x10 Wipes per case
Art. no: 5568

Year End Price! 1 190 SEK
List price 1 640 SEK

AGMA

All offers, except the Promega offer, are valid until December 31, 2024, and can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please refer to promo code "YearEnd24" in order to get the right price when ordering.



Are frozen cell lines slowing down your research? You're not alone!

Many of our partners in cell line research, diagnostics, and drug development face a common challenge: the complex, time-consuming thawing process of frozen cells.

This process can:

- Take days or even weeks
- Require specialized expertise
- Be prone to errors
- Create unexpected delays
- Lead to misunderstandings between suppliers and clients

The solution? "Ready-to-use" cells and cell lines (including organoids!)

By eliminating the thawing process and sending the samples in the portable CO2 incubator, you can:

- Increase productivity
- Improve cost efficiency
- Enable precise time management
- Ensure consistent quality from day one

Our customers working with cells, organoids, 3D cell structures, and tissue patches have seen transformative results. Not only does this approach optimize the entire workflow, but it also ensures satisfaction on both sides.

Please contact us for more information.

CELLBOX™
LIVE CELL SHIPMENT

NanoPhotometer maintenance and calibration free

UV-Vis Spectrophotometers made in Germany, totally maintenance and calibration free and can be operated via a mobile app.

- Scanning in the area from 200 - 900 nm can be done in 2.5 seconds
- The built-in vortex makes sure the sample is homogenized properly
- Autosampling

Quantify and qualify your nucleic acid, protein, bacterial culture and many other liquid samples with the most convenient micro volume spectrophotometer ever built.



N120

Measure 12 samples (1-2 µl) in parallel. Read time: 1.7 s/sample, 20 s/run for 12 samples.

N120-TOUCH
Art no: N120
Year End Price! 178 000 SEK



NP80

For measuring 0.3-2 µl sample or from a cuvette - all-in-one instrument. Read time: 3.5 s/sample

NP80-TOUCH
Art no: NP80-TOUCH
Year End Price! 120 000 SEK



N60

For measuring small samples (0.3-2 µl).

NP60-TOUCH
Art no: N60-TOUCH
Year End Price! 104 000 SEK

Welcome the new Implen family member OD600!






- Modern, stand-alone, battery powered device with direct operation via the touchscreen display
- Designed with simplicity in mind, the OD600 features an intuitive interface for easy operation
- Fast and easy determination of cell density at 600 nm
- Measure OD600 values up to 4 A

IMPLEN OD600
Art no: IMPLENOD600-EU
Year End Price! 13 000 SEK
List price 16 500 SEK

The Device is Pre-programmed with Intuitive OD600 and McFarland Applications



All offers, except the Promega offer, are valid until December 31, 2024, and can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please refer to promo code "YearEnd24" in order to get the right price when ordering.





FastPrep-24™ 5G Instrument

Get fast, bench-top sample homogenization with the convenience of 72 pre-defined protocols and space for 12 user-defined protocols. Compatible with 2 mL, 4.5 mL, 15 mL and 50 mL tubes. Offers simultaneous homogenization of up to 24 samples within 40 seconds.



FastPrep-24™ Classic Instrument

Affordable and easy to use, the FastPrep-24 Classic enables consistent lysis for reproducible extraction of DNA, RNA, proteins. Offers simultaneous homogenization of up to 24 samples within 40 seconds



FastPrep-96™ Instrument

Maximize productivity with this high-throughput bead beating system that delivers on performance, speed, and reproducibility. Compatible with 2 mL, 4.5 mL, 15 mL, 50 mL tubes, 250 mL flasks and 96 tubes in a plate.



SuperFastPrep-2™

A hand-held sample disruption instrument designed to quickly and efficiently lyse biological samples at the bench or in the field. Compatible with 2 mL tubes.



FastPrep Lysing Matrix tubes

Lysing Matrices to use with the FastPrep come in 16 unique matrix compositions to ensure reliable homogenization of any sample type preparation from day-to-day and sample-to-sample. Validated worldwide with 3,000+ Lysing Matrix specific publications.



FastPrep Adapters

Designed for sample sizes ranging from 2 mL to 250 mL, our adapters increase the flexibility of your sample preparation.

Small ULT Freezers when space is limited

The freezers use eco-friendly, HCFC-free Nature R refrigerant.



Benchtop freezer (-80°C)

WxDxH: 630 x 675 x 650 mm
Volume: 37 L

Art no: N123000
Year End price! 58 830 SEK
List price: 78 441 SEK



Slim upright freezer (-80°C)

WxDxH: 630 x 370 x 1420 mm
Volume: 256 L

Art no: N123002
Right now! 79 700 SEK
List price: 113 861 SEK



Comfort Rack **Standard Rack:**



Blizzard Ultralow Temperature Freezers utilize non-toxic hydro-carbon (HC) refrigerants that are eco-friendly alternatives to CFC, HCFC, and HFC fluorocarbons, which are linked to ozone layer damage.

ULT freezer for many samples (-80°C)

WxDxH: 1145 x 980 x 1980 mm
Volume: Upright 828 L
Capacity: 600 cryo boxes

Art no: NU-99828JE-N
Year End Price! 132 000 SEK
List price: 199 566 SEK

Racks fitting the 828L freezer

Offer valid when buying 24 racks
(which corresponds to a full freezer)

Comfort Rack
Sliding drawers, reduced operating weight
Rack fits 25 cryoboxes, height 50 mm
278x683x140mm

Art no: TE24254
Year End Price! 1 870 SEK/each
(x24 rack: 44 880 SEK)
List price: 2 772 SEK /st (x24 rack: 66 528 SEK)

Standard Rack:
Standard rack, fits 25 cryoboxes, height 50 mm
278x683x139mm

Art no: TE23154
Year End Price! 1 190 SEK/each
(x24 rack: 28 560 SEK)
List price: 1 738 SEK/st (x24 st rack 41 712 SEK)

NORDIC LAB NUAIRE TENAK

New supplier of freezers and climate cabinets

Apollo Laboratory Technology have their production in the Netherlands and are specialized in thermal equipment; machines used to heat, cool or condition. Besides freezers and climate chambers they also offer refrigerators, pass through refrigerators, ovens and ultra low freezers.

The freezers have an isolation made of 70 mm (the fridge/freezers have an isolation of 60 mm) ecological polyurethane foam in the walls making them a sustainable, energy-efficient solution but also ensures optimal maintenance of the desired temperature of your samples. The freezers come as standard with both audible and visual alarms when the door is open or the temperature is too high or low. Optional power failure alarm can also be selected.



Freezer -10…-30°C
Capacity: 80 L
2 shelves
Dimensions: 650x645x853 mm
Noise level: 51 db

Art no: SLTM 100
Year End Price! 30 600 SEK
List price: 38 200 SEK



Freezer -10…-30°C
Capacity: 260 L
4 shelves
Dimensions: 540x760x1994 mm
Noise level: 53 db

Art no: SLTM 300
Year End Price! 32 600 SEK
List price: 40 700 SEK



Fridge/Freezer +2..+10°C/-10..-30°C
Capacity: 120 L+120 L
4 shelves
Dimensions: 520x751x1994 mm
Noise level: 53 db

Art no: ALC SLTM 300
Year End Price! 43 000 SEK
List price: 55 300 SEK



Climate cabinet 700

The 700 climate cabinets have a temperature range of 0°C to + 70°C and an air humidity range of 10% to 98%. The cabinets are very precise in their control, accurate in achieving the values and very uniform.


Capacity: 580 L
4 (adjustable) shelves
Dimensions: 740x870(824)x1994 mm
Noise level: 55 dBa

Art no: CLH 700
Year End Price! 140 000 SEK

APOLLO
Laboratory Technology

Cryogenic ULT -150°C freezers available

The cryogenic -150°C chest freezer from PHC provide the most uniform storage temperatures for long-term cryopreservation solutions. The freezers are available in 128 L (art. MDF-1156-PE) and 231 L (art. MDF-C2156VAN-PE).



- Optimum sample storage due to use of patented space saving VIP Plus panels within the freezer cabinet
- Compressors are specifically designed for ultra-low temperature applications
- Specially designed LCD control panel provides a clear view of the temperature and gives notification in case of abnormalities



Charlotte Otéma
070 889 17 59
charlotte.otema@nordicbiolabs.se



Feven Solomon
073 373 21 09
feven.solomon@nordicbiolabs.se



Jens Götharson
073 668 29 90
jens.gotharson@nordicbiolabs.se



Peggy Roth
073 668 29 91
peggy.roth@nordicbiolabs.se

All offers, except the Promega offer are valid until December 31, 2024, and can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please refer to promo code "YearEnd24" in order to get the right price when ordering. Nordic Biolabs cannot be held responsible for any printing errors and/or price alterations.

nordic biolabs
www.nordicbiolabs.se

Nordic Biolabs AB, Box 7293, 187 14 Täby,
Tel vxl 08-630 85 00, Fax 08-630 09 05, info@nordicbiolabs.se