NORDIC BIOLABS AB PRESENTS

Products for Life Science

PROMO CODE: YEAREND26-e

4/2026

NEW Thermo-Shaker

That heats and cools – now with upgraded platforms **and** optional heated lid to minimize condensation during temperature-controlled workflows.

TS-100C Pro, Thermo-Shaker with cooling

Art. No.: BS-010178-A01

Right Now!
22 500 SEK*

**Blocks and lid are sold separately*



nordic biolabs

www.nordicbiolabs.se info@nordicbiolabs.se 08-630 85 00

NEW

Sartorius Light Pack™ – The Responsible Tip Box

*Designed to deliver a measurable reduction in environmental impact, the **Sartorius Light Pack™** combines sustainability with laboratory convenience. Compared to conventional single-tray plastic racks, the plastic content has been reduced by **74%**, helping laboratories minimize plastic waste.*

The Light Pack features **FSC-certified cardboard**, a **recyclable plastic tray**, and a window made from **cellulose acetate**, supporting responsible material choices throughout the packaging design.


Optilid
Hinged lid which can be opened and closed with one hand

Foldable
At end of use, reducing waste volume by up to 50%

Cardboard box with Plastic Tray
Plastic reduced by 74%




Cellulose Acetate Window
to see the number of tips in the box



Robust and Rigid Design
For stability during tip picking

Safety Corner and Seamless Bottom
Cleanliness and reduction in contamination risk



After use, the box folds flat for space-efficient storage, easy stacking, and seamless integration into existing recycling streams.

Awarded with the **My Green Lab® ACT® Ecolabel.**

my green lab. ACT Ecolabel

reddot winner 2026

SARTORIUS

Offer on Pipettes and Tips

YEAR-END OFFER

Buy a **20, 100, or 200µl Tacta** or **120µl Picus 2** pipette and get **1 pack of Optifit Light Pack** pipette tips sterile or unsterile for FREE.

Art. No.	Description
LH-729030	Tacta Mechanical Pipette, 1-ch, 2-20 µl
LH-729050	Tacta Mechanical Pipette, 1-ch, 10-100 µl
LH-729060	Tacta Mechanical Pipette, 1-ch, 20-200 µl
LH-747041	Picus 2 electronic pipette, 1-ch, 5-120 µl
LH-747341	Picus 2 electronic pipette, 8-ch, 5-120 µl
LH-747441	Picus 2 electronic pipette, 12-ch, 5-120 µl



Choose 1 pack from below for free:

Art. No.	Description
LH-G790200	Optifit 200µl, Light Pack (10x96)
LH-G790201	Optifit 200µl Pre-st, Light Pack (10x96)

To claim your offer, add promo code "YEAREND26-e" to your order when purchasing in the webshop.

Invest in Confidence. Invest in *Sartorius*.



Exclusive offer for laboratories seeking accuracy, reliability, and compliance.

In modern laboratories, every milligram counts. Whether you work in pharmaceutical research, biotechnology, diagnostics, academia, or quality control, the **Sartorius Entris™ II Essential Analytical Balance** delivers the precision, stability, and reliability required for confident results every day.

Art. No.	Description	Capacity / Readability	List Price (SEK)	Price Right Now! (SEK)
BCE64I-1S	Analytical balance	60 g / 0.1 mg+ isoCAL	25 670	20 298
BCE224I-1S	Entris II Essential balance	220 g / 0.1 mg + isoCAL	28 490	22 533

Top Features

- LED touch technology for intuitive navigation and use
- Powerful and precise monolithic weigh cell technology exclusively available with Sartorius balances
- Automatic calibration and adjustment documentation with isoCAL
- 12 integrated applications expand usability and application potential
- Modern connectivity options for simple data transfer and analysis

Popular Applications

- Routine weighing tasks during research, production, processing, and more
- High measurement precision and accuracy during routine weighing tasks in laboratories, research centers, production and processing facilities, and more

Contact us today for a quotation, demonstration, or application support.

SARTORIUS

All offers are valid until December 31, 2026 (if not stated otherwise), and can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please add promo code "YEAREND26-e" to your order to receive the right price.



Small Instruments for a Small Price!

Art. No.	Description	List Price (SEK)	Price Right Now! (SEK)
BS-010152-AAG 	**MR-1 Mini Rocker/Shaker** Small shaker for regulated gentle rocking motion ideal for mini gel destaining after electrophoresis, conducting Northern, Southern, and Western blot analysis.	8 210	5 410
BS-010203-AAG 	**V-1 Personal vortex** Two modes: 1. Continuous operation; 2. Impulse operation (activated by pressing the cap with the tube's bottom).	3 880	2 550
BS-010151-AAG 	**3D, Sunflower Mini-Shaker** Mini-shaker is designed for mixing blood samples, for minigel staining and destaining, sample washing, and blot hybridization reactions.	6 310	4 155
BS-010128-AAI 	**PST-60HL-4 Plate shaker-thermostat** Thermo–Shaker designed for shaking and thermostating 4 standard 96–well microplates.	26 890	17 730
BS-010202-AAA 	**FVL-2400N Mini-centrifuge/Vortex Combispin** For simultaneous mixing and separation of samples, using centrifuge and mixing modules, located on the common spin–module. Rotation speed: 3500 rpm, max RCF 700 x g.	6 890	4 360
BS-040108-AAG 	**FTA-1, Flask-trap aspirator, volume 1L** Designed for aspiration/removal of residual liquids (buffers, alcohols, etc.) from cell culture vessels, test tubes or microtubes, in cell culture or molecular biology (e.g. DNA/RNA precipitation).	7 220	4 670

Fast Single Tube Reader

Fast small benchtop scanner for reading 1D coded racks and tubes or 2D coded tubes. Simply plug into a computer and get the data straight into Excel without the need for any software installation.



Micronic Tube Reader DT520

Art. No.: MP55129

Right Now! 11 200 SEK



MICRONIC
INNOVATION IN SAMPLE STORAGE



Get Excited – Fluorescence 96!

Compact, single-channel fluorometer brings big-league performance right to your benchtop, adding much-needed capacity for your assays without taking over the workbench.

- **Practically unlimited dynamic range:** No need to worry about saturation!
- **Top or bottom reads:** Flip the device setup effortlessly
- **Fast, parallelized reads:** 8-channel detection gets you from plate to results faster
- **Open Format:** Your familiar fluorescence assays in standard 96-well format
- **Small:** Dimensions: 9.7 x 16.0 x 5.8 cm (W x L x H), Weight: 723 g

Fluoresence 96

Right Now! 83 000 SEK

byonoy

All offers are valid until December 31, 2026 (if not stated otherwise), and can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please add promo code "YEAREND26-e" to your order to receive the right price.



Something Cool for Your Budget

Promega's *popular promotion is back — buy 4 identical products and get 25% discount!*

Winter is here. So are the savings. Let your budget breathe easy this winter and **save 25% on Promega reagents.**

But do not wait too long. The campaign runs from **November 1st** until **December 15th.**

This promotion is valid on all Promega's standard catalog products, including Maxwell kits, FuGENE products, instruments, and Glo-products. Please note that it cannot be combined with other promotions or special agreements.

How to receive your free product?

- **Step 1:** Order at least 4 of the same product
- **Step 2:** Put promo code "**4f3**" on your order.
- **Step 3:** Please ensure your order is placed before **Tuesday 15 December 2026.**
- **Step 4:** We will send you 4 products for the price of 3, saving you 25% on your year-end budget. Plus, you can easily repeat the process if you need more products!

Just make sure your order includes promo code 4f3, and the more you order, the more you save — while supplies last.*

A selection of Promega products:




PCR Related Products
for amplifying DNA and RNA sequences for applications in scientific research, medical diagnostics, and forensic analysis.



RNasin® Ribonuclease Inhibitors
for RNA isolation, RT-PCR, in vitro transcription or any situation where RNase contamination is a concern.



Automated Nucleic Acid Extraction
using the Maxwell chemistry.



Cell-Based Assays
for measuring cell viability, stress, apoptosis, and death.

** For questions, special quotes, or larger quantities, contact Charlotte Otéma at charlotte.otema@nordicbiolabs.se.*

RNAsure
– Protect RNA from the Moment of Collection



RNAsure® *is a specialized buffer system designed to help stabilize RNA at room temperature.*

Samples collected in RNAsure® can be stored/transported as follows without compromising RNA quality:

- 37°C: up to 1 day
- 25°C: up to 1 week
- 4°C: up to 1 month
- -20°C to -80°C: long-term storage

Art. No.	Description	Price Right Now! (SEK)
2010100100	RNAsure, 100 mL	1 560
2010100500	RNAsure, 500 mL	5 850

INVITEK diagnostics

Analyze Open Chromatin Regions at Genome-Wide Scale

ATAC-Seq *is a rapid assay that allows analysis of epigenetic profiles across the genome by identification of regions that have open or accessible chromatin states.*

Because of the assay's speed, simplicity, and applicability to a wide range of sample types, ATAC-Seq has become a commonly-used epigenetic assay, and it can serve as a gateway to further, more detailed epigenetic analyses.

In the ATAC-Seq assay, intact nuclei from cell or tissue samples are treated with a hyperactive Tn5 transposase mutant. This transposase is able to simultaneously tag the target DNA with sequencing adapters and fragment the DNA in a process termed tagmentation.

ATAC-Seq Kit

Art. No.: 53150-AF

Right Now! 20 % Discount



ACTIVE MOTIF®
Tools to Analyze Nuclear Function

All offers are valid until December 31, 2026 (if not stated otherwise), and can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please add promo code "YEAREND26-e" to your order to receive the right price.

DISCOVER OUR PCR CABINETS






Premium Product Class

Premium Product Class

	UVC/T-M-AR	UVT-B-AR	UVC/T-AR	UVT-S-AR (double size)
Walls material	Rear: stainless steel Sides and front: glass	Rear: stainless steel Sides: steel with chemical resistant powder coating Front: glass	Plexiglas: Polymethyl methacrylate Altuglas EX	Rear: stainless steel Sides and front: glass
Working Area	645 x 190 mm	645 x 190 mm	645 x 165 mm	1210 x 190 mm
Outer Dimensions	700 x 580 x 555 mm	700 x 580 x 555 mm	720 x 535 x 555 mm	1250 x 600 x 590 mm
Inlet for Power Cords	Inlet for power cords & 1 built-in socket	Inlet for power cords & 1 built-in socket	Inlet for power cords	Inlet for power cords & 3 built-in socket
Art. No.	BS-040104-A30	BS-040109-A20	BS-040102-A13	BS-040107-A07
Year End Price! (SEK)	**30 000**	**25 600**	**23 800**	**41 700**
List Price (SEK)	40 060	35 120	32 680	57 330

Premium Product Class *UVC/T-M-AR and UVT-S-AR are defined as Premium Class offering one additional year of warranty (2+1)*

For work in clean, controlled, and contamination-free environments when working with nucleic acids and conducting PCR workflows.

- **High-Powered UV Sterilization** — Effectively destroys contaminants and prevents cross-contamination
- **Continuous Protection with UV Recirculators** — Ensures ongoing sterilization during PCR setup
- **Digital Display with Programmable Timer** — Allows users to track UV lamp life and set sterilization schedules
- **UV Safety Features** — Equipped with UV protection film to ensure safe operation while using high-intensity UV sterilization



Eliminate DNA & RNA from Surface Prior PCR Reaction Preparation

PDS-250, ready-to-go formulation in a spray bottle, 250ml

Art. No.: BS-040107-DK

Right Now!
460 SEK/bottle


biosan

Still Waiting on LC-MS to QC Your mRNA?

CELLSCRIPT™ ***EZ-QC™ mRNA Assay Kits*** *bring quantitative mRNA QC to your own benchtop — no outsourcing, specialized equipment, or multi-day wait.*

Using standard polyacrylamide gel electrophoresis (PAGE), EZ-QC™ assays let you evaluate 5′-capping efficiency, Cap 1 efficiency and 3′-poly(A) tail length — validated directly against LC-MS:

- Capping and Cap 1 efficiency results within ±5% of LC-MS
- Poly(A) tail length measurements closely matching LC-MS
- Same-day results, right at your own benchtop
- Lower sample input — as little as 5 picomoles, vs. 150-250 picomoles typically needed for LC-MS



EZ-QC™ mRNA Cap 1 Efficiency Assay Kit

Quantitates the percentage of Cap 1 vs Cap 0 caps of synthesized mRNA

Art. No.: ONE240910

Right Now! 5 485 SEK

CELLSCRIPT™
RNA for Translation in Cells

Special Price on *Implen* OD600®!



Implen OD600®


Art. No.: IMPLENOD600-EU


Right Now! 15 100 SEK

IMPLEN

- Modern, stand-alone, battery powered device with direct operation via the touchscreen display
- Designed with simplicity in mind, the OD600 features an intuitive interface for easy operation
- Fast and easy determination of cell density at 600 nm
- Measure OD600 values up to 4 A

The Device is Pre-programmed with Intuitive OD600 and McFarland Applications





All offers are valid until December 31, 2026 (if not stated otherwise), and can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please add promo code "YEAREND26-e" to your order to receive the right price.

Fridges, Freezers, and Racks

Pharmacy Fridges

Range from 36 L to 543 L

A next generation pharmacy refrigerator featuring the latest IntelliCold® 24-Range technology with touch screen controller. Available with either solid or glass door and key lock or digital lock.

Searching for another model? Contact us for a quote.



LABCOLD

Art. No.	Description	Price Right Now! (SEK)
RLDF0124	Pharmacy Fridge 36L, Solid Door, Lock	14 785
RLDF0224	Pharmacy Fridge 66L, Solid Door, Lock	16 299
RLDF0524	Pharmacy Fridge 150L, Solid Door, Lock	20 800
RLDF1524	Pharmacy Fridge 440L, Solid Door, Lock	41 185

Chest or Upright ULT Freezers

Range from 37 L to 828 L

Searching for another model? Contact us for a quote.



Small Under Counter Freezer

Upright 91 L, -86°C Freezer ULT U100
WxDxH: 595 x 585 x 825 mm

Art. No.: N123001

Year End Price!
57 300 SEK

List Price: 90 030 SEK



Large Capacity Upright Freezer

828 L, -86°C freezer
WxDxH: 1145 x 980 x 1980 mm
Capacity: 600 cryoboxes

Art. No.: NU-99828JE-N

Year End Price!
138 100 SEK

List Price: 199 566 SEK



We Provide Racks for All Freezers on the Market

Contact us with your specific need and we will help you to find the right rack.

NORDIC LAB NUAIRE TENAK



Label with Confidence.

Work with Certainty.

*Whether you are managing clinical samples, research materials, cryogenic storage, or quality control processes, the **M210-LAB, M510** and **M611** provide the durability and functionality required for demanding laboratory environments.*


Art. No.	Description	Recommendation for Laboratory Environments	Price Right Now! (SEK)
311319 	**M210-LAB Portable Label Printer**	Best for simple, durable laboratory labeling without software.	1 816
322723 	**M510 Label Printer EU + Brady Workstation Suite**	Best balance of laboratory use, flexibility, and PC integration.	6 070
322920 	**M611 Label Printer EU + Workstation Suite**	Best choice when wireless connectivity, mobility, and the widest range of label materials are required.	8 561

Feature	M210-LAB	M510	M611
Print Speed	10 mm per sec	38.1 mm per sec	50.8 mm per sec
Cutter type	Manual cutter	Manual cutter	Automatic cutter
Label Options	90+ label options	250+ standard parts	800+ standard parts
Label Colors	Monocolour	Monocolour	Monocolour
Connectivity options	None	USB	Bluetooth, Wi Fi, USB
PC Compatibility	No	Yes	Yes
Custom Labels Available	No	No	Yes
Interface type	Keyboard only	Full QWERTY Keyboard, PC	Smartphone/PC controlled

Contact us today to learn more, request pricing, or schedule a hands-on product demonstration.

BRADY

All offers are valid until December 31, 2026 (if not stated otherwise), and can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please add promo code "YEAREND26-e" to your order to receive the right price.



Time to Replace FBS with a Synthetic Product?

*We are thrilled to offer **SEREXIS™**, a product designed to replace FBS, addressing the growing concerns about sustainability and animal welfare. By choosing a synthetic product you also eliminate the variability seen in a biological product like FBS.*

Use this as a replacement, so no change of protocols! We suggest you try first to verify that it works in your cell culture experiments.

If You Prefer Sticking to FBS – Test Batch Before Buying!

Test and reserve – we offer testing of 2-3 different batches when reserving a minimum of 20 bottles of each batch. For larger needs we offer storage of the bottles once you have decided.

Contact Charlotte Otéma for you sample.

cytiva CAPRICORN SCIENTIFIC



Stock Up on D-PBS

Now for an Even Better Price!

Dulbeccos PBS (1x) w/o Ca, Mg, Phenol Red (500 mL)

Art. No.: PBS-1A

Right Now! 47 SEK/bottle

The campaign price is valid when buying one box = 20 bottles.

Component	Concentration mg/L
KCl	200
KH2PO4	200
NaCl	8000
Na2HPO4 (anhydrous)	1150

CAPRICORN SCIENTIFIC



Incusafe Multigas Incubators:
The Reliable Solution for Labs

*The **Incusafe multigas incubators** are equipped with a unique Zirconia O_2 sensor and an IR2 CO_2 sensor, which ensure fast and accurate measurements of oxygen and CO_2 concentrations.*

InCuSaFe®, a copper-enriched stainless steel interior, inhibits microbial growth, reduces maintenance requirements, and lowers the total cost of ownership.

Art. No.	Description	Price Right Now! (SEK)
MCO-50M-PE	Multigas incubator, 50 L, Dual CO_2 IR + Zirconia O_2 Sensor	129 450
MCO-50MUV-PE	Multigas incubator, 50 L, Dual CO_2 IR + Zirconia O_2 Sensor, UV kit	137 970

NEW

Cellbox™ GO

***Cellbox™ Go Ground**, designed for short- and mid-distance transport of live cells, tissues, organoids and other cell-based materials, keeps temperature and carbon dioxide (CO_2) levels precisely controlled while your samples are on the move.*




Cellbox™ Go Ground
Art. No.: 10000540


New Core Edition
Affordable, Functional, and Easy to Use

The New Core Edition offers the essential functionalities of the Cellbox™ Go Ground with a more affordable configuration. Temperature and CO_2 settings are pre-set to your specifications before delivery, so it is ready to provide controlled conditions during cell hand-offs between labs.

How It Works:



01 – Select your temperature



02 – Select your CO_2 level




03 – Pre-configured for you



04 – Ready for your transport

All offers are valid until December 31, 2026 (if not stated otherwise), and can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please add promo code "YEAREND26-e" to your order to receive the right price.



TPP SlideMax™ Chamber Slides

TPP **SlideMax™** enables the direct cultivation of cells on a standard **Cyclo-Olefin Polymer (COP) microscope slide**, making it an ideal solution for immunofluorescence, staining applications and advanced microscopy workflows where exceptionally low autofluorescence is essential.

Art. No.	Description	CS	List Price (SEK)	Price Right Now! (SEK)
71011	Slidemax 1 chamber, lid, PS/COP, 93x33x32 mm	5	1 050	540
71012	Slidemax 2 chambers, lid, PS/COP, 93x33x32 mm	5	1 050	540
71014	Slidemax 4 chambers, lid, PS/COP, 93x33x32 mm	5	1 050	540
71016	Slidemax 6 chambers, lid, PS/COP, 93x33x32 mm	5	1 050	540
71018	Slidemax 8 chambers, lid, PS/COP, 93x33x32 mm	5	1 050	540

Nordic Biolabs Suppliers

Innovative epigenetics tools and services for research and drug discovery

Centrifuges for GMP and GLP environments

Disinfection and cleaning for lab and cleanrooms

Equipment for heating, cooling, and airconditioning in laboratory and industrial test environments

Large system for cryo-/biobanks

Ice-free sample preparation and freezing

Antisera, media, bacteria definition

Lab instruments with small footprint

Labels, printers, and spill protection

Aborbance readers with small footprint

Cell culture media and serum

Portable CO_2 incubators for live cell transport

New generation bioreactors

RNA for translation in cells

Swabs, sample collection, and transport

Cell media and serum for research and drug discovery

Sterilization control

Tubes, sample containers, and storage

Equipment and associated products for microbiology and cell culture applications

Biochemicals, bio-reagents, and antibodies

Tubes, racks, reagent reservoir, centrifuges, magnetic stirrer etc.

Wipes for cleanroom applications

Nanophotometers and cuvette photometers

Instruments and IVD labeled DNA/RNA kits

Ethanols and desinfection

Laboratory, pharmacy, vaccine refrigeration and blood product storage

Library prep kits for RNA

Laboratory glassware such as microscope slides, cover glasses, and counting chambers

Microbiological media and supplements

Tubes in 96 well format for storage

Biochemicals, lysing matrix, homogenizers for DNA, RNA and proteins from any tissue

Low temperature lab freezers

Osmometers, automated cell culture analyzers, and systems for liquid measurement

CO_2 incubators, safety cabinets, and ULT-freezers

Chromatography consumables, sample prep for bioanalysis, and proteomics

Nitrogen gas generator and hydrogen gas generator systems for LC-MS, GC and GC-MS

Freezers, autoclaves, and incubators

Environmental monitoring

Cell- and molecular biology reagents, multimode readers, and extraction systems

Cleanroom and industrial safety products

Accurate, automated cell count and analysis

Pipette tips, filtration and purification, pipettes, balances, centrifuges, water systems, and air samplers

Nitrile gloves for personal protection

Sterilization control

Cuvettes and certified reference material

Consumables and accessories for HPLC/UHPLC, GC, and SPE

Racks for all freezer brands and models

Sterile cell culture plastics, Pipette boy (Turbofix)

Contamination control and environmental monitoring

Sales & Portfolio Managers



Charlotte Otéma

+46 8 630 85 03
charlotte.otema@nordicbiolabs.se



Peggy Roth

+46 8 630 85 24
peggy.roth@nordicbiolabs.se



Feven Solomon

+46 8 630 85 16
feven.solomon@nordicbiolabs.se



Marit Hoppe

+46 8 630 85 22
marit.hoppe@nordicbiolabs.se

Nordic Biolabs AB is environmentally certified according to ISO 14001. In this mailing, we have made a conscious decision to use a thinner paper that we believe will be better for the environment.

Nordic Biolabs AB
Ritarslingan 16A
187 66 Täby
Sverige

Phone +46 8 630 85 00
Email info@nordicbiolabs.se
Website www.nordicbiolabs.se / www.nordicbiolabs.dk

nordic biolabs

All offers are valid until December 31, 2026 (if not stated otherwise), and can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please add promo code "YEAREND26-e" to your order to receive the right price. Nordic Biolabs cannot be held responsible for any printing errors and/or price alterations.