Pipette Doctor at Nordic Biolabs

We service and calibrate all brands and models.



We now calibrate
multichannel pipettes at
our laboratory in Täby!

nordic biolabs

www.nordicbiolabs.se pipett@nordicbiolabs.se +46 8 630 85 05

About the Pipette Doctor

The Pipette Doctor offers professional calibration, servicing, and repair of pipettes — all to ensure the quality of your work and extend the lifespan of your pipettes.

The Pipette Doctor is part of Nordic Biolabs, a privately owned distributor of laboratory products, instruments, and technical services that has been active on the Swedish market since 1970.

The Pipette Doctor — Established in 2018
In 2018, we opened our pipette calibration laboratory, Pipette Doctor™, at our headquarters in Täby, just north of Stockholm. This made it possible for Swedish customers to have their pipettes calibrated locally.

Since 2021, the laboratory has been accredited according to ISO 17025, ensuring that our calibrations meet international standards for quality and competence.

Quality and Sustainability at the Core

Our certifications and accreditations include:

- **ISO 9001** – Quality
- **ISO 14001** – Environment
- **ISO 17025** – Pipette Calibration

Meet the Team!



Lisa Tibbelin
Manager Pipette Doctor
+46 8 630 85 26
lisa.tibbelin@nordicbiolabs.se



Albin Söderblom
Service Technician



Joanna Tagai
Service Technician



Mikael Selin
Service Technician



Nicoleta Ursu
Service Technician



Welcome to Our New Colleague!

The Pipette Doctor continues to grow, and we are excited to welcome our new service technician — Svitlana Semeniuk — to the team!

NEW SERVICE

We Now Calibrate Your Multichannel Pipettes in Täby!

During 2025, we have expanded our laboratory with the Speedcal Mobile System. This means you can now have your multichannel pipettes calibrated locally at our lab in Täby.



With our new multichannel balance system, we can now offer fast calibration of multichannel pipettes directly on-site in Täby. For you as a customer, this means shorter lead times, less transport handling, and a greater security.



> **Please note:** Calibration of multichannel pipettes up to 10 µl, as well as accredited calibration of multichannel pipettes, is performed by our partner Sartorius in Finland.

We welcome both new and existing customers to take advantage of this improved service.

Pipette Calibration & Service

We provide a full range of services in calibration, maintenance, repair, and training — all to help you maintain the quality of your work and extend the lifespan of your pipettes.

Our process:

1. **Choose your service level***
2. **Send in your pipettes** *(we also offer pipette pickup services)*
   Clean your pipettes with 70% ethanol and complete the decontamination form. It can be downloaded from our website — www.nordicbiolabs.se/pipettdoktorn

   Pack your pipettes, the completed decontamination form, and tips (if desired), and send them to:

   Nordic Biolabs AB
   Att. Pipettdoktorn
   Ritarslingan 16A
   187 66 TÄBY
3. **Maintenance**
   Your pipettes are cleaned, lubricated, and reassembled according to the manufacturer's instructions.
4. **Repairs**
   If needed, we perform service and repairs on your pipettes — this is included in the price (excluding spare parts). A cost estimate will be provided if the cost of spare parts exceeds SEK 600 per pipette.
5. **Calibration**
   Performed according to current MPE tolerances in ISO 8655:2022, using either original tips or customer-supplied tips (6 tips per channel). For multichannel pipettes, all channels are tested.
6. **Report**
   After service is completed, you will receive your pipettes along with a service and calibration report.

***Service levels**

Level	Includes	Pre-test	Maintenance	Measurements	Report	ISO 17025 Certificate
1	Precalibration	✅	⛔	5 or 10 measurements on 2 or 3 volumes	✅	⛔
2	Maintenance & Calibration (5.2)	⛔	✅	5 measurements on 2 volumes	✅	⛔
3	Maintenance & Calibration (10.2)	⛔	✅	10 measurements on 2 volumes	✅	⛔
4	Maintenance & Calibration (10.3)	⛔	✅	10 measurements on 3 volumes	✅	⛔

Accredited Pipette Calibration & Service

We offer accredited calibration in accordance with ISO 17025, ensuring that our measurement results are traceable, reliable, and meet international quality standards.

Why choose accredited calibration?

☑ Recommended for organizations with high quality standards, certifications, or regulatory requirements.
☑ Performed according to ISO 17025 with traceability to international standards.
☑ Includes complete documentation and certificate.

	Pre-test	Maintenance	Measurements	Report	ISO 17025 Certificate
Accredited Calibration (ISO 17025)	✅	✅	10 measurements on 3 volumes	✅	✅

The Pipette Academy — Training & Seminars

Would you like to improve the precision, accuracy, and comfort of your pipetting technique? Our seminars combine theory and hands-on training to give you the tools you need for reliable results every time.

During the seminar, you will learn:


- Pipetting techniques
- Ensure accuracy and precision at every step
- Identifying and preventing common errors
- Maintenance and calibration according to ISO 8655
- Ergonomics and risk reduction in daily work



Feel free to contact us
— we are happy to tell you more!

Nordic Biolabs — We Are More Than Our Service

Nordic Biolabs offers a wide range of laboratory products and is a leading provider in areas such as transfusion, cell culture, microbiological sampling for diagnostics, environmental monitoring in production, and molecular biology. Discover our selection of laboratory products and instruments in areas such as...



Life Science



Transfusion



Diagnostics



Cleanroom

... to name a few. Explore our full range on our website!

Contact Us

If you have any questions regarding the Pipette Doctor, pricing, service agreements, or pipette pickup, feel free to contact us.

Pipette Doctor

Address	Ritarslingan 16A, 187 66 Täby
Opening Hours	08.00-16.30 (08.00-16.00 summer time)
Website	www.nordicbiolabs.se/pipettdoktorn
Tel Direct	+46 8 630 85 05
Email	pipett@nordicbiolabs.se

nordic biolabs

Nordic Biolabs AB, Ritarslingan 16A, 187 66 Täby
Tel vxl +46 8-630 85 00, info@nordicbiolabs.se

Calibrates and maintains
all makes – all models – all volumes