NORDIC BIOLABS AB PRESENTS

Products for Life Science

Promo code: NBLQ3


3/2024


Introducing the Implen OD600®



- Modern, stand-alone, battery powered device with direct operation via the touchscreen display
- Designed with simplicity in mind, the OD600 features an intuitive interface for easy operation
- Fast and easy determination of cell density at 600 nm
- Measure OD600 values up to 4 A

The Device is Pre-programmed with Intuitive OD600 and McFarland Applications





nordic biolabs

www.nordicbiolabs.se info@nordicbiolabs.se ordertel 08-630 85 10



New supplier

WE WELCOME APOLLO LABTECH!

Providing thermal equipment used to heat, cool or condition

Apollo Labtech also supply and develop equipment for cleaning, (quick) drying or humidification among others. These can be standard products as well as innovative custom solutions. Working closely with their partners, they strive every day to exceed their customers` expectations and continue to amaze them.



Apollo´s climate cabinets

- Wide range in capacity - 300 to 1400 liters
- Temperature range +0°C - +70°C
- Humidity range 10% - 98%
- Moisture fluctuation ± 4%

Product range

- Refrigerators
- Pass through refrigerators
- Refrigerators with suction
- Ovens
- Freezers
- Ultra low freezers
- Climate cabinets
- Climate cabinets with CO_2
- Climate rooms
- Syn extraction

APOLLO Laboratory Technology

All offers are valid until October 31, 2024, and can not be combined with other offers or agreements with Nordic Biolabs AB.
All prices are ex VAT. Please refer to promo code "NBLQ3" in order to get the right price when ordering.




NEW CENTRIFUGE

Small size centrifuge available in a ventilated or refrigerated version

Choose from 11 different rotors for centrifuging for example 48 microtubes, 6 x 50 mL tubes or 6 microplates.

No need to touch the centrifuge to start your cycle. A simple gesture and you can start, stop and open your lid.

- Unique illuminated display, visible at 7 meters
- A green circle or blinking red circle confirm whether or not the cycle is running smoothly
- **Change rotor without tools** - All rotors are equipped with the exclusive Clickspin™ system, allowing you to quickly and easily change the rotor



Angular Rotor
24x2mL - 45° - Sealed
AFI-LINA24-2



Angular Rotor
6x50mL - 35°
AFI-LINA6-50



Angular rotor PCR
8x8ml - 45° - Sealed
AFI-LINAPCR8-8

Description	
AFI-CI030V-E	LINA Centrifuge Ventilated (V)
AFI-CI030R-E	LINA Centrifuge Refrigerated (R)
Specifications	Angular Rotor
Max Speed	23 000 rpm / 39 625 g
Max loading	6 x 50 ml
Dimensions	28x31x40 cm (V)/28x35x64 cm (R)
Weight	21 (V) /40 (R) kg

ACT label consumables at Nordic Biolabs with lower environmental impact

The ACT Ecolabel program by My Green Lab aims to address the need of both scientists and procurement specialists for clear, third-party verified information about the environmental impact of laboratory products. The ACT Ecolabel, which stands for Accountability, Consistency, and Transparency (ACT), is a standardized measure of a product's environmental impact, assisting laboratories in making informed, sustainably conscious choices when comparing products from different manufacturers. The lower the score, the lower the impact on the environment in terms of manufacturing, energy and water use, packaging, and end-of-life

ACT. The Environmental Impact Factor Label — EU

Sartorius Optifit Pipette Tips Refill Tower 350µl
Kajaani, Kainuu, Finland
SKU 790352

Environmental Impact Scale
1 Decreasing Environmental Impact 10

Manufacturing	
Manufacturing Impact Reduction	1.0
Renewable Energy Use	Yes
Responsible Chemical Management	1.0
Shipping Impact	1.0
Product Content	5.0
Packaging Content	1.0
User Impact	
Energy Consumption (kWh/day)	N/A
Water Consumption (liters/day)	N/A
Product Lifetime	6.8
End of Life	
Packaging	1.0
Product	7.9
Innovation	
Innovative Practices	-1.0
Environmental Impact Factor:	**23.7**
Label Valid Through:	March 2026

act.mygreenlab.org



Refill Pipette Tips

The Optifit Tip Refill Tower offers a variety of tip volumes in a space-saving package for a wide range of laboratory applications. Refill your empty Single Trays from the Optifit Tip Refill Tower to reduce plastic consumption and logistics emissions. The Refill Tower packaging has earned the ACT label.

Refill packs provide an effective method for reducing plastic waste in laboratory, with the added benefit of being offered in pre-sterilized variants to meet the stringent sterility requirements of the more purity demanding applications.

ACT. The Environmental Impact Factor Label SARTORIUS

1.4 mL Storage Tubes Internal Thead

The Micronic 1.40 mL tubes were independently verified by My Green Lab, and now carry the ACT label.

The Micronic 1.40ml storage tubes with internal thread and flat outer bottom are provided with a code surface that cannot be separated from the transparent tube due the unique 2K injection molding technique. As a result, tube codes can never wear or fall off. The 1.40ml tubes are standard available in the 96-well format Micronic 96-4 Rack.

Combine with the recycled Micronic EFC 96-3 Rack for even better choice for the environment.





MICRONIC INNOVATION IN SAMPLE STORAGE

ACT. The Environmental Impact Factor Label — EU

Micronic 1.40ml 2D Data-Matrix coded Screw Cap tubes V-bottom precapped with Grey Screw Cap in 96-4 Rack with high cover (Barcoded A1-H1 side)
Lelystad, Netherlands
SKU MP52551-Z20

Environmental Impact Scale
1 Decreasing Environmental Impact 10

Manufacturing	
Manufacturing Impact Reduction	1.0
Renewable Energy Use	No
Responsible Chemical Management	1.0
Shipping Impact	1.0
Product Content	10.0
Packaging Content	1.0
User Impact	
Energy Consumption (kWh/day)	N/A
Water Consumption (liters/day)	N/A
Product Lifetime	6.0
End of Life	
Packaging	4.0
Product	8.8
Environmental Impact Factor:	**32.8**
Label Valid Through:	June 2023

act.mygreenlab.org

All offers are valid until October 31, 2024, and can not be combined with other offers or agreements with Nordic Biolabs AB.
All prices are ex VAT. Please refer to promo code "NBLQ3" in order to get the right price when ordering.



PEAK Scientific - our newest supplier of chromatography related products

Most chromatography users are familiar with PEAK Scientific as a leading provider of high-performance nitrogen gas generator and hydrogen gas generator systems for LC-MS, GC and GC-MS.

We will work closely together with Peak Scientific to secure our customers the best gas generator solution.
Please contact us for any questions and requests regarding PEAK Scientific or other of our chromatography products.

Contact
Iben Kofoed Piper
+45 9384 5256
iben.kofoedpiper@nordicbiolabs.dk

Move chemical bottles in a safer way

Safety Bottle Carriers designed to hold various-sized bottles in a safe manner and protect the bottle from breakage. The bottle carriers are made of flexible plastic that is chemical resistant and autoclavable.



The carriers are available in green & blue and include an ergonomic handle to prevent hand fatigue.
The handle swings out of the way for easy placement of the bottle in the carrier.

- Features a recessed marking area at the top of the carrier
- Includes a centering ring in the base
- Dishwasher safe
- Autoclavable at 121°C
- Comes in 2 different sizes (1.5 liters and 4.5 liters)



NEW!

The Quintix Pro Balance from Sartorius

With the Quintix® Pro, Sartorius replaces the Secura® and Quintix® product lines. The new balance Quintix® Pro is designed to deliver unparalleled performance with premium features that enhance ease of use and adaptability in different laboratory environments.*

(*) The Secura® and Quintix® lines will be discontinued by the end of 2024, but Quintix® semi-micro balances will remain in the portfolio.

- Easy access to menu option with a large, high-resolution full-touch graphic display with modern upgraded user interface
- Internal calibration
- Quick access to the integrated digital manual
- Seamless connectivity features include interface options, easy data transfer, and GLP | GMP compliant printout
- **Certification for use in cleanrooms grade A and B, and vacuum, argon, and nitrogen environments**
- Chemically-resistant surfaces made of PBT, glass and stainless steel support routine cleaning protocols

Reduced environmental impact from the manufacturing plant to the bench due to;

- Climate-conscience production in addition to environment-friendly packaging
- Power-saving operation with stand-by or switch-off models
- Digital user manual minimizes paper consumption, while allowing quick access to information

SARTORIUS

All offers are valid until October 31, 2024, and can not be combined with other offers or agreements with Nordic Biolabs AB.
All prices are ex VAT. Please refer to promo code "NBLQ3" in order to get the right price when ordering.

Weighing tables from VISU

VISU's weighing tables are designed for equipment which is sensitive to vibration, such as scales. The tables are of large mass and have damping elements made of rubber which eliminate disturbances effectively.



- Easy to clean.
- Available at standard and customized heights
- Available with electrically adjustable height

Description	Product code (numbers indicate standard dimensions width-depth-height [mm])	Height [mm]	Weight approx. [kg]
Weighing table with electrically adjustable height	EVD-1200-550-S	Adjustment range from 650 to 940	170
Weighing table with electrically adjustable height	EVD-900-550-S	Adjustment range from 680 to 970	160
Weighing table with electrically adjustable height	EVD-550-550-S	Adjustment range from 680 to 970	100
Weighing table with a metal base	EVA-550-550-750	750	100
Weighing table with stone base	EVB-900-550-750	750	300
Weighing stone on a worktop	EVC-300-300-30	30	10





Weigh boats in different shapes and material

Weigh boats in **polystyrene** come in flat-bottom square, hexagon and diamond shapes.

Aluminium boats have crimped edges for sturdiness and includes a finger tab.

SmartBoats offers a green solution to sample weighing made of FSC certified **paper** and when disposed in the standard waste stream they will naturally biodegrade (Eco-OUT).



High Speed Magnetic Stirrers

Compact stirrers with stainless steel working surface. Speed control 0-3000 rpm.



MS-3000

Made of strong steel and painted with powder enamel, which is resistant to acids and alkali. Supplied with a cylinder-shape magnetic stirring bar (6 × 25 mm) encapsulated in PTFE for universal use.



MMS-3000

Equipped with a detachable stand for supporting varoius sensor elements (temperature, pH and others) inside the stirred liquid. Strong magnets hold the driven magnetic element firmly in the magnetic clutch. Stirring is performed without undesirable heating and noise.

biosan
Medical-Biological Research & Technologies

	MS-3000	MMS-3000
Article number	BS-010301-AAF	BS-010305-AAF
Stirring volume up to (H_2O)	5 L	20 L
Max length of magnetic stirring bar	50 mm	70 mm
Working plate size	110 x 100 mm	Ø160 mm
Price right now!	**2 500 SEK**	**2 800 SEK**
List price	3 580 SEK	3 950 SEK

NEW KIT!

Purify plasmid DNA in less than 15 min

The kit employs an advanced silica-based membrane technology, presented in the form of a user-friendly spin column.

- Recover up to 20 µg of high-quality plasmid DNA, making it suitable for various molecular biology applications such as cloning and in vitro transcription
- The InviSorb® streamlined procedure eliminates the need for additional extraction or precipitation steps, providing high-quality plasmid DNA that is immediately ready for use in downstream applications



INVITEK diagnostics

Art no	Description	Right now! SEK
1010140300	InviSorb® Spin Plasmid Miniprep Kit, 250 preps	3 300

All offers are valid until October 31, 2024, and can not be combined with other offers or agreements with Nordic Biolabs AB.
All prices are ex VAT. Please refer to promo code "NBLQ3" in order to get the right price when ordering.

Calibrate your pipettes regularly!

Depending on how much your pipettes are used and of how many people are using them this can vary from every 3rd month to once a year, but it will ensure that the pipettes perform at highest precision and accuracy over time.

By planning the maintenance and calibration in good time before problems arise you avoid potential errors in future work. Make sure you regularly clean and take care of the pipettes also between the time of calibration.
If you note that you have an inconsistency in for example the volume delivery you likely need to service and calibrate your pipette.


Pipette Doctor Täby service all makes, models and volumes

We offer:

- The full package which means that service and repair also are included in the price (ex spare parts)
- Pick-up of pipettes (the regions around Stockholm, Uppsala, Gothenburg, Linköping and Malmö / Lund)
- Service, calibration and repair of multichannel pipettes is done at Sartorius in Helsinki
  (also applies to accredited calibration according to ISO 17025)

If you have questions about the Pipette Doctor service, please contact us.

pipett@nordicbiolabs.se
08 630 85 05





FBS Minis

Contamination? *Aliquoting FBS into smaller units is a common cause of contamination. Controlled clean room conditions eliminate this concern. Additionally, potential sources of contamination, such as pipettes and the risk of "reflux drops" during aliquoting, are non-existent with FBS Minis.*

Enhanced reproducibility for your research - Reproducibility remains a challenge in cell culture experiments. Due to the strictly controlled manufacturing process, every parameter is meticulously inspected and of course you can test your batch the batch before deciding.

Eco friendly & resource-efficient With FBS Minis, you contribute to the reduction of plastic waste, as there is no need for additional pipettes or tubes for aliquoting (saving of consumables approx. 30 %). In addition, you save working time, thereby freeing up important resources to be used elsewhere!

CAPRICORN SCIENTIFIC

Art no	Description	Pack size
10-FBS-12F	FBS Minis, South American	10 x 50 ml

Searching for Affordable Serum Alternatives?

Engineered HyClone Serum from Cytiva - Reliable results, better availability, better economy! Developed for better performance and specialized use.


- **General alternative to FBS**
  HyClone bovine growth Serum, US origin, an affordable and viable alternative to Fetal Bovine Serum
- **For Hybridoma growth**
  HyClone Fetal Clone I serum: less variation, low e-tox and Hb-levels
- **For CHO cells and derivatives**
  HyClone Fetal Clone II serum: improved growth, less variation
- **For fibroblasts and many other cells**
  HyClone Fetal Clone III serum, US origin: better growth, less variation
- **For antibody production in hybridomas**
  HyClone Alpha Calf Fraction media supplement, a cost-effective medium supplement



FREE SAMPLE and lot reservation!
charlotte.otema@nordicbiolabs.se

cytiva

All offers are valid until October 31, 2024, and can not be combined with other offers or agreements with Nordic Biolabs AB.
All prices are ex VAT. Please refer to promo code "NBLQ3" in order to get the right price when ordering.



Trypsin Digestion in as Little as 60 Minutes

The Rapid-Digestion Trypsin and Rapid-Digestion Trypsin/Lys-C Kits are designed to shorten protein digestion times to 60 minutes versus the typical 4–18 hours required for trypsin digestion.

Both kits contain three components: i) protease (Trypsin or Trypsin/Lys-C Mix); ii) protease Resuspension Buffer; and iii) Rapid-Digestion Buffer optimized for faster protein digestions.

Trypsin digestion with these kits follows a simple-to-use protocol that is both fast and efficient. The protocol is flexible, can accommodate a large range of sample volumes and protein concentrations and requires no special laboratory equipment or off-line desalting. The entire sample preparation procedure is performed in as little as 60 minutes.

Promega

Art No	Description
VA1060	Rapid Digestion Kit---Trypsin, 100 µg
VA1061	Rapid Digestion Kit---Trypsin/Lys-C, 100 µg

Lysyl Endopeptidase

Lysyl Endopeptidase®, Mass Spectrometry Grade (Lys-C)
Proteins can be identified by mass spectrometry analysis of the peptides produced by in-gel digestion, and further information regarding post-translational modifications can be obtained.

Lysyl Endopeptidase®, Mass Spectrometry Grade (Lys-C) is a freeze dried product that retained sufficient activity for in-gel digestion and packed in very small quantities for convenience purposes.


FUJIFILM Value from Innovation Wako

Art no	Description	Listprice SEK	Price per piece when buying at least 2 packs, SEK
125-05061	Lysyl Endopeptidase Mass spectrometry grade, 20 ug x 5/fp	18 200	14 560



Freezer dimensions
Outer: 595 x 630 x 825 mm
Inner: 430 x 426/205 x 625/395 mm

Stackable low temperature freezers



By using a coupler, the freezers can be stacked safely on top of each other.

When space is limited and/or sample volumes at first are small you can choose a flexible system which can be expanded over time.

The freezers are designed as undercounter units but are ideal for stacking.
The small 91 L freezer, is available in two different temperatures and can be stacked in pairs of two, in different combinations using a coupler.

NORDIC LAB

Art no	Description	List price SEK	Right now! SEK
N122001	XLT U100 Upright Freezer 91L -°65	65 500	52 400
N123001	ULT U100 Upright Freezer 91L -86°C	83 160	66 500



Charlotte Otéma
070 889 17 59
charlotte.otema@nordicbiolabs.se



Feven Solomon
073 373 21 09
feven.solomon@nordicbiolabs.se



Jens Götharson
073 668 29 90
jens.gotharson@nordicbiolabs.se



Peggy Roth
073 668 29 91
peggy.roth@nordicbiolabs.se

All offers are valid until October 31, 2024, and can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please refer to promo code "NBLQ3" in order to get the right price when ordering. Nordic Biolabs cannot be held responsible for any printing errors and/or price alterations.

nordic biolabs
www.nordicbiolabs.se

Nordic Biolabs AB, Box 7293, 187 14 Täby,
Tel vxl 08-630 85 00, Fax 08-630 09 05, info@nordicbiolabs.se