Centrifuges from AFI

A selection of centrifuges in different sizes, ventilated, refrigerated or heated.

- Up to 18 interchangeable rotors are available
- Click-spin rotors – easy change of rotor without using a tool
  New feature: A green light confirms that the rotor is successfully locked
- New feature: Heated version available - From ambient to 37°C in only 9 minutes



Click Spin with visual lock indicator



LIA Microcentrifuge

Vent.	Refrig.	Heated	Max Speed Swing out rotor	Max Speed Angular rotor	Max loading (Swing-out/Angular)	Nb of Rotors
x	x		N/A	15000 rpm / 21380g	NA /18x5 ml	8 rotors



LINA Microcentrifuge

Vent.	Refrig.	Heated	Max Speed Swing out rotor	Max Speed Angular rotor	Max loading (Swing-out/Angular)	Nb of Rotors
x	x		N/A	23000 rpm / 39625g	NA/ 6x50 ml	11 rotors



SIRENA (0,6 L)

Vent.	Refrig.	Heated	Max Speed Swing out rotor	Max Speed Angular rotor	Max loading (Swing-out/Angular)	Nb of Rotors
	x		4700 rpm / 3725g	18000 rpm / 31594g	4x150 ml / 8x100 ml	17 rotors



ISA (0,6 L)

Vent.	Refrig.	Heated	Max Speed Swing out rotor	Max Speed Angular rotor	Max loading (Swing-out/Angular)	Nb of Rotors
x	x	x	4500 rpm / 3396g	6700 rpm / 6331g	4x150 ml / 6x85 ml	11 rotors



LISA (2,5 L)

Vent.	Refrig.	Heated	Max Speed Swing out rotor	Max Speed Angular rotor	Max loading (Swing-out/Angular)	Nb of Rotors
x	x	x	5100 rpm / 5816g	18000 rpm / 31514g	4x625 ml / 8x100 ml	18 rotors



LOREENA (4 L)

Vent.	Refrig.	Heated	Max Speed Swing out rotor	Max Speed Angular rotor	Max loading (Swing-out/Angular)	Nb of Rotors
x	x	x	4600 rpm / 4731 g	16000 rpm / 28335g	4x1000 ml / 6x250 ml	15 rotors

Version 2, February 22, 2024

nordic biolabs

www.nordicbiolabs.se info@nordicbiolabs.se ordertel 08-630 85 10